

easyProphy®

Air Polising & Cleaning System



Before



After



Ergonomic Design

The easyProphy is designed to give you the maximum comfort with well balanced handpiece.



Flexible

easyProphy turns 360 degrees. The handpiece connection is specially designed to rotate freely. The handpiece can be repeatedly autoclaved at up to 135°C.



Swivel and Quick Release

Simple pull back the quick release ring to release the easyProphy from the coupling.



Cleaning with taste

With special mint and lemon flavoured cleaning powder. Dental treatment tastes even better to your patients.



Working Principle

The compressed air flows into the powder container when the foot control is pressed. It is mixed uniformly with the powder, and the mixture flows to the tip. At the same time, water is supplied separately from air and powder.

There are three water ports located at the tip to create water sheath. It surrounds the powder to prevent dusty operative field.

Intended usages

- Removing of plaque and stains.
- Polishing after calculus removal.
- Cleaning of enamel before etching for subsequent fissure sealing.
- Preparing of cavities for improved bonding between enamel and filling materials.
- Cleaning of implant surfaces.

easyProphy leads to a new era of tooth cleaning and polishing. It operates in a very simple way that directly attaches to the handpiece couple or tubing. It fits into 4,5,6 hole tubing.

The anti-clogging design will maximize your efficiency and reduce your frustration.

Code	Model	Description
EP3001	easyProphy 200/M4	Standard 4 hole connection
EP3003	easyProphy 200/QD	Quick Disconnect, Swivel
EP3009	H1	easyProphy Handpiece with Nozzle
EP3010	N1	easyProphy Nozzle
PR1506	O-ring	O-ring Set (5Pcs/Pack)
PV1100	Powder	Mint flavor (130g/btl)
PV1101	Powder	Lemon flavor (130g/btl)