Result Reference Guide

Unprocessed

Accept Regular Cycles

Accept Extended Cycles

Reject

NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

Explanation of Symbols

Product designed for use STEAM with Steam sterilization cycles.

Batch number.

CE mark.

Expiration Date.

Authorized representative in the European Community.

Manufacturer.

Intended Use Table

Model	Autoclave / Steam Cycles		
	Gravity Displacement		
CI5-P1	121 °C	132 °C	135 °C
	30 minutes	15 minutes 25 minutes	10 minutes
	Dynamic Air Removal (Vacuum Assist)		
	121 °C	132 °C	135 °C
	NA	4 minutes	3 minutes
	Minimum Stated Values (SV) for Chemical Integrators		
	121 °C	132 °C	135 °C
	16.5 minutes	2.0 minutes	1.2 minutes

EN Chemical Integrator

for Steam Sterilization

Indications for use

United States

The integrator Beyes® Masteri® CI5-E1 is designed to chemically react over time with the critical parameters of steam sterilization cycle within a specified tolerance. The integrating indicator strip is intended to be placed in each pack, pouch, container, tray or other containment device to function as an independent monitor of critical parameters for the following sterilization cycles: Gravity Displacement: 121°C for 30 minutes, 132°C for 15 minutes, 132°C for 25 minutes, 135°C for 10 minutes; Dynamic Air Removal (Vacuum Assist): 132°C for 4 minutes, 135°C for 3 minutes. SV121°C/16.5 min., SV132°C/2.0 min, SV135°C/1.2 min. See Intended Use Table for reference.

Outside the United States

Masteri® CI5-E1 moving front chemical integrator is designed for monitoring Steam sterilization processes between 118°C and 138°C.

Device description

Masteri® CI5-E1 Chemical Integrator (Type 5 according to ISO 11140 -1:2014 standard) are single-use moving front chemical indicators for monitoring Steam sterilization processes. These indicators ensure an adequate control of the effectiveness of steam sterilization processes by monitoring all critical parameters. Masteri® CI5-E1 integrator indicator consists of a special paper wick and a steam and temperature sensitive chemical pellet contained in a paper /film /foil laminate. The chemical pellet melts and migrates as a dark bar along the paper wick. The migration is visible through a zone marked ACCEPT or REJECT, thus indicating whether sterilization conditions were met. The extent of migration depends on steam quality, time and temperature. The ACCEPT result is reached when a theoretical spore population reaches its kill time, indicating integration condition has been reached. This condition is calibrated with the kill time of a 106 Geobacillus stearothermophilus spore population, calculated in BIER (Biological Indicator Evaluator Resistometer).

Precautions

Do not expose Masteri® CI5-E1 to EO, Dry Heat, or any sterilization process other than Steam. Do not reuse. Do not reuse the sterilizer until the integrator reaches integrating condition.

Instructions for use

- Carefully open the package and remove only the number of indicators needed. Reseal the package tightly.
- Place a Masteri® CI5-E1 indicator in each pack, peel pouch, container system or tray to be steam sterilized in the most inaccessible area to sterilizing agent.
- Perform de sterilization cycle.
- After the sterilization process has finished, the dark bar should have reached the ACCEPT zone. If the dark bar has not reached the ACCEPT zone, a REJECT result is indicated and the items in the pack, peel pouch, container system, or tray were not exposed to proper steam sterilization conditions. These items should be reprocessed according to established procedures and a new CI5-E1 has to be used.

NOTE: In some extended or wet steam cycles, the dark bar can display discoloration in some zones of the reading window, imposing a difficulty for result's interpretation (see Result Reference Guide for reference). Despite this, regarding cycle integration parameters, this process can be considered as accepted in case of extended cycle. However, a damp appearance of the CI5-E1 indicator after the sterilization process is evidence of a wet cycle. In this case, the cycle should be rejected, even if

the dark bar reached the ACCEPTED zone, and Beyes® recommends planning an equipment's maintenance.

NOTE: If any serious incident occurs in relation to the device, it should be reported to Beyes® and the competent authority of the State in which the user is established.

Monitoring frequency

Follow facility Policies and Procedures which should specify a chemical indicator monitoring frequency compliant with professional association recommended practices and/or national guidelines and standards. As a good practice and to provide optimal patient safety, Beyes® recommends that every sterilization load be monitored with an appropriate chemical indicator.

Storage

Indicators should be stored in a dry place, protected from light at temperatures between 10-30°C, 30-80% relative humidity. Do not wet. Do not store near sterilizing agents.

Shelf life

Chemical Indicators have an expiration date of 5 years from the date of manufacture when stored in recommended conditions. Both the expiration date and batch number of the product are specified on the package label. Do not use indicators after their expiration date.

Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal

Discard chemical indicators after use along with paper waste, according to your country's healthcare and safety regulations.

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