
## Unique Point Chemical Integrator for Steam Sterilization



### Result Reference Guide

#### Unprocessed


#### Failures

One or more parameters of sterilization were not achieved. The materials cannot be considered sterile and must be processed again. These results are usually obtained when steam quality is poor or adequate vacuum has not been generated.

#### Processed and Satisfactory

**NOTE:** Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

#### Explanation of Symbols

 Product designed for use STEAM with Steam sterilization cycles.

 Batch number.

 CE mark.

 Expiration Date.

 Authorized representative in the European Community.

 Manufacturer.

#### Intended Use Table

Model	Autoclave / Steam Cycles		
	Gravity Displacement		
CI5-P1	121 °C	132 °C	135 °C
	30 minutes	15 minutes 25 minutes	10 minutes
	Dynamic Air Removal (Vacuum Assist)		
	121 °C	132 °C	135 °C
	NA	4 minutes	3 minutes
	Minimum Stated Values (SV) for Chemical Integrators		
	121 °C	132 °C	135 °C
	16.5 minutes	2.0 minutes	1.2 minutes

## EN Unique Point Chemical Integrator

for Steam Sterilization

#### Indications for use

##### United States

The integrator Beyes® Masteri® CI5-P1 is designed to chemically react over time with the critical parameters of Steam sterilization cycle within a specified tolerance. The integrating indicator strip is intended to be placed in each pack, pouch, container, tray or other containment device to function as an independent monitor of critical parameters for the following sterilization cycles: Gravity Displacement: 121 °C for 30 minutes, 132 °C for 15 minutes, 132 °C for 25 minutes, 135 °C for 10 minutes; Dynamic Air Removal (Vacuum Assist): 132 °C for 4 minutes, 135 °C for 3 minutes. SV121 °C/16.5 min., SV132 °C/2.0 min, SV135 °C/1.2 min. See Intended Use Table for reference.

##### Outside the United States

Masteri® CI5-P1 unique point chemical integrator indicator is designed to monitor Steam sterilization processes between 121 °C and 135 °C.

#### Device description

Masteri® CI5-P1 Chemical Integrator (Type 5 according to ISO 11140-1:2014 standard) are single-use unique point chemical indicators that consist of paper strips printed with indicator ink. They have been designed for Steam sterilization processes, ensuring an adequate control of the effectiveness of sterilization processes. Yellow ink was developed to turn to black when a theoretical spore population reaches its kill time, indicating integration condition has been reached. This condition is calibrated with the kill time of a 10<sup>6</sup> Geobacillus stearothermophilus spore population, calculated in BIER (Biological Indicator Evaluator Resistometer).

Indicators 100 % free of toxic heavy metals.

#### Precautions

Do not expose Masteri® CI5-P1 to EO, Dry Heat, or any sterilization process other than Steam. Do not reuse. Do not reuse the sterilizer until the integrator color reaches integrating condition.

#### Instructions for use

- Place a chemical integrator Masteri® CI5-P1 within each package, bag or tray that will be sterilized by Steam. In rigid containers place an CI5-P1 in every corner, or at least two diagonally opposite corners, to ensure that steam has reached the most inaccessible areas.
- Perform the sterilization cycle.
- After the sterilization process has finished, remove the indicators from the sterilizer and analyze the results. For all the chemical indicators exposed to Steam, the ink must turn to the black reference color when the integration conditions are reached. For reference color, please refer to Result Reference Guide.

**NOTE:** if any serious incident occurs in relation to the device, it should be reported to Beyes® and the competent authority of the State in which the user is established.

## Monitoring frequency

Follow facility policies and procedures which should specify a chemical indicator monitoring frequency compliant with professional association recommended practices and/or national guidelines and standards. As a good practice and to provide optimal patient safety, Beyes® recommends that every sterilization load be monitored with an appropriate chemical indicator.

## Storage

Indicators should be stored in a dry place, protected from light at temperatures between 10-30 °C, 30-80 % relative humidity. Do not wet. Do not store near sterilizing agents.

## Shelf life

Chemical Indicators have an expiration date of 5 years from the date of manufacture when stored in recommended conditions. Both the expiration date and batch number of the product are specified on the package label and in each unit of the product. Do not use indicators after their expiration date.

Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

## Disposal

Discard chemical indicators after use along with paper waste, according to your country's healthcare and safety regulations.

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Beyes Dental Canada Inc.
Lotus NL B.V.