

AirLight

High Speed Handpiece Collection



-
- Built-in Generator
 - Direct-LED

Innovation

1. Built-in Generator

A micro generator is built into the handpiece to power its LED. It eliminates the need for a fiber optic system.

2. Direct-LED

A macro angle LED is optimally positioned on the head to provide wider and brighter illumination.

Solutions Giver

Say Goodbye to....

- Fragile Light Bulbs
- Expensive Fiber Optic Systems
- Yellowish and Dimming Light



AirLight



Regular Fiber Optic

Money Saver \$\$\$

AirLight help doctors to
save...

\$1000+ for installing fiber optic systems

\$100+ for replace broken optic tubings

\$100+ for light bulbs





Key Features

AirLight

Triple Water Spray

Effectively cools the entire surface of the bur that avoid overheating and provides better patient comfort.

Ceramic Ball Bearing

It is specially designed for smooth and quiet operation. It also provide superior torque and long lifespan.

Micro-Tex Coating

The latest surface enhancement technology creates a micro-textured and scratch proof surface that lowers the possibility of slippery during operation.



What doctors say...

AirLight



"It was very light and easy to maneuver".

Frank Kuzmin, DMD CT, USA



"The AirLight...seems like an efficient cutting instrument acting as a natural extension of my hand."

Frank M. Glushefski, DMD PA, USA



"AirLight is the solution for better visibility and no more fiber optic system break down. The torque and balance are impressive "

Allen Gaw, DDS TX, USA

Beyes 

Proven to be the leader!!



We are proud that AirLight is in the TOP 100 list. An annual highlight in Dentistry Today, will take a look at the TOP 100 products during the past year based on lead generation. It focuses on some of the newest and most innovative products in the dental profession



AirLight receive 4 stars from THE DENTAL ADVISOR after 1 year clinical trial. Since 1984, THE DENTAL ADVISOR has provided dental professionals worldwide with concise, accurate, and objective information using in-vivo and in-vitro results on the performance of products and equipment based on clinical and laboratory research.



AirLight received "Best Evaluated Product" in 2012. There were 25 dentists involved in evaluating this product for half year term in USA. AirLight received a spectacular result, which 96% evaluators would recommend this products to their colleagues



The "Recommended Product" awards products that have received high score in their independent product review process. Dental Product Shopper's independent evaluators review hundreds of products each year. Our AirLight qualified as a "Recommended Product"

Technical Specification

- Standard and Mini Head
- Push Button Auto Chuck
- Ceramic Ball Bearings
- Mini Torque :17 Watts
- Standard Torque : 21 Watts
- Noise Level: 59 decibels
- Intensity: 25,000 Lux
- Autoclavable: 135°C
- 350,000-400,000rpm
- 3-ports Coolant Spray
- Swivel and Quick Disconnect

Order Information



Beyes QD

HP3014, M800-S/QD, Standard Head

HP3015, M800-M/QD, Mini Head



Kavo* Connection

HP3035, M800-S/K, Standard Head

HP3036, M800-M/K, Mini Head



NSK * Connection

HP3042, M800-S/N, Standard Head

HP3043, M800-M/N,, Mini Head



Star * Connection

HP3053, M800-S/ST, Standard Head

HP3053, M800-S/ST, Mini Head



W& H * Connection

HP3021, M800-S/W, Standard Head

HP3022, M800-M/W, Mini Head



www.beyes.ca

info@beyes.ca

order@beyes.ca

Beyes Dental Canada Inc
9-595 Middlefield Road,
Toronto, Ontario, M1V 3S2, Canada
1-416-293-2868 1-416-800-7475

AirLight