

# CANAVIEW<sup>®</sup>

## INTRAORAL CAMERA

An excellent communication tool that improves understanding and interaction with patients. It provides visual evidence that you need to educate patients and make more accurate diagnosis. Boasting a true autofocus and intuitive camera lighting system, the dental digital photography camera requires no manual adjustments and is easily shared between different examination rooms and chairs.

### New Features

### 3 Zoom Options • 3 Light Intensity Options

*(Only for Canaview V2)*

- True Auto-Focus
- Anti-Fog Design
- Responsive Capture Button
- 5 MP High Resolution Images
- Intuitive Camera Lighting System
- No Manual Adjustments Required
- Direct USB Connection to Computer
- Super Light-Weight



### Complimentary Image Management Software:

Abdent, Adstra, Apixia, Apteryx, CaptureLink, Clio Soft, Cuve Hero, DBSWIN, Dentrix Ascend, Dentrix Image, Dexis, Diagora, Eaglesoft, EVAsoft, ImageFX, Kodak, MiPACS, Mogo, Open Dental, Romexis, Schick, Sidexis, Sopro Imaging, Tiger view, Umbie, VaTech EzDent, VisiQuick, Visix, VixWin.....



## Key Specifications

<b>Focus Type</b>	Auto
<b>Focus Range</b>	4mm - 60mm
<b>Connection Type</b>	USB to computer
<b>Image Quality</b>	5 million pixels
<b>Resolution</b>	1600 x 1200
<b>Camera Length</b>	20cm
<b>Weight</b>	35g
<b>Cable Length</b>	2m
<b>Power</b>	DC 5V+/-0.5V(1.5A)
<b>Operating System</b>	Windows 10/8/7/NXP/Vista (64/32 bit)
<b>Warranty</b>	12 months



## Order Information

- Code:** EV2060  
**Model:** Canaview Intraoral Camera
- Code:** EV2080  
**Model:** Canaview V2 Intraoral Camera
- Code:** EVP008  
**Model:** Canaview Holder



**Beyes Dental Canada Inc.**  
 23-595 Middlefield Rd  
 Scarborough, ON M1V 3S2  
 Canada  
**Tel: 1-855-603-1888**  
**Email: info@beyes.ca**  
**Website: www.beyes.ca**



**Lotus NL B.V.**  
 Koningin Julianaplein 10, Ie Verd,  
 2595AA, The Hague,  
 Netherlands

### Your Dealer for Beyes Products



Scan the QR code  
 for more information

Printed in Canada  
 File #: ENB057  
 Rev. 1/12.04.24

© Beyes Canada Inc 2024. All Rights Reserved.  
 Specifications are subject to change without notice.  
 All brands are holders of their respective trademarks.  
 All Beyes products are manufactured in strict compliance  
 with Health Canada regulations ensuring quality.  
 Manufacturing facility is ISO 13485 certified under the  
 Medical Device Single Audit Program (MDSAP).