

Maxso[®]

High Speed Handpiece

The Maxso series high speed handpiece is designed to meet expectations in modern dental procedures. High quality components and comfortable grip with Micro-tex finish provide superior handling that you desire. Furthermore, the well aimed triple water spray will provide excellent coolant effects. Maxso series comes with variety of connections; as a result, you simply attach the Maxso to your existing system.



Triple Water Spray

Effectively cools the entire surface of the bur that avoids overheating and provides better patient comfort.





Micro-Tex Coating

The latest surface coating technology creates a micro-textured and scratch proof to enhance handling and durability.



Compatibility

Compatible with a great variety of coupling systems and adapts great to every existing treatment units.





Technical Data

Standard Head

Mini Head









Head Size(Width x Height)	12.5mm x 13.1mm	10.8 x 13.0mm
Drive Air Pressure(PSI)	40-45 PSI	40-45 PSI
Power(WATT)	21	19
Speed(RPM)	350,000	400,000
Rotary Instruments	FG burs according to EN ISO 1797-1:1995	FG burs according to EN ISO 1797-1:1995
Ø / max. length(MM)	1.6 / 25	1.6 / 25
Chuck Type	Push Button	Push Button
Water Spray	Triple	Triple
Midwest 4 Hole Connection		
 	CODE HP2008 MODEL M200-S/M4	CODE HP2009 MODEL M200-M/M4
Borden 2 Hole Connection		
 	CODE HP2012 MODEL M200-S/B2	CODE HP2013 MODEL M200-M/B2
Beyes® PD Connection		
 	CODE HP6022 MODEL M200-S/PD	CODE HP6012 MODEL M200-M/PD
NSK® Connection		
 	CODE HP2036 MODEL M200-S/N	CODE HP2037 MODEL M200-M/N



Technical Data

Standard Head

Mini Head

Head Size(Width x Height)	12.5mm x 13.1mm	10.8 x 13.0mm
Drive Air Pressure(PSI)	40-45 PSI	40-45 PSI
Power(WATT)	21	19
Speed(RPM)	350,000	400,000
Rotary Instruments	FG burs according to EN ISO 1797-1:1995	FG burs according to EN ISO 1797-1:1995
Ø / max. length(MM)	1.6 / 25	1.6 / 25
Chuck Type	Wrench Chuck	Wrench Chuck
Water Spray	Single	Single
Midwest 4 Hole Connection		
 	CODE HP2001 MODEL SA20-S/M4	CODE HP2000 MODEL SA20-M/M4
Borden 2 Hole Connection		
 	CODE HP2004 MODEL SA20-S/B2	CODE HP2003 MODEL SA20-M/B2
Chuck Type	Push Button	Push Button
Water Spray	Single	Single
Midwest 4 Hole Connection		
 	CODE HP2106 MODEL M200E-S/M4	CODE HP2105 MODEL M200E-M/M4
Borden 2 Hole Connection		
 	CODE HP2109 MODEL M200E-S/B2	CODE HP2108 MODEL M200E-M/B2