

# SLIMAX-C Plus<sup>®</sup>

## LED Curing Light System



The newest Slimax-C Plus Curing Light System provides an innovative built-in digital radiometer that gives you an instant output measurement. You will be more confident than ever when calculating the curing time. This system is also equipped with an LED display panel to give you the control that you need. It is powered by a rechargeable Samsung Lithium-ion battery that provides up to 850 cures per charge. It also generates powerful output up to 2000 mW/cm<sup>2</sup> ± 10%. The Slimax C-Plus comes with a one year warranty.

- Built-in Radiometer
- Instant Output Measurement

### Durable

Solid aluminum body sets it apart from all its plastic competitors

### Portable

Cordless design with a long-lasting battery

### Powerful

Generates powerful output of 2000mW/cm<sup>2</sup>(±10%)

### Beautiful

3 colors give you a décor choice with a personal touch



## Key Specifications

<b>Power Supply</b>	Lithium-ion battery (3.7V 2200mAh)
<b>Wavelength Range</b>	430-480 nm
<b>Wavelength Peak</b>	455 nm ± 10 nm
<b>Light Intensity</b>	2000 mW/cm <sup>2</sup> ± 10% (independent of battery power level)
<b>Total Exposure Time (fully charged battery)</b>	850 times / 10-second cures
<b>Charging Time (depleted battery)</b>	Approximately 3 hours
<b>Dimensions</b>	Diameter 23 mm (0.92"), Length 254mm (10")
<b>Weight</b>	160g (0.35 lb) (includes light guide)

## Order Information



Code: **SM1003P-B**  
Model: **Slimax-C Plus Blue**



Code: **SM1003P-PI**  
Model: **Slimax-C Plus Pink**



Code: **SM1003P-S**  
Model: **Slimax-C Plus Silver**

Code: **DB1024**  
Model: **SLMS-1**  
Description: Full Curing Light Sleeve, 500PCS/Box

**Beyes**  
CANADA

**Beyes Dental Canada Inc.**  
23-595 Middlefield Rd  
Scarborough, ON M1V 3S2  
Canada  
Tel: 1-855-603-1888  
Email: info@beyes.ca  
Website: www.beyes.ca



**Lotus NL B.V.**  
Koningin Julianaplein 10, Ie Verd,  
2595AA, The Hague,  
Netherlands

Your Dealer for Beyes Products



Scan the QR code  
for more information

Printed in Canada  
File #: ENB006R4  
Rev.2/09.16.24

© Beyes Canada Inc 2024. All Rights Reserved.  
Specifications are subject to change without notice.  
All brands are holders of their respective trademarks.  
All Beyes products are manufactured in strict compliance  
with Health Canada regulations ensuring quality.  
Manufacturing facility is ISO 13485 certified under the  
Medical Device Single Audit Program (MDSAP).