

# CareMaster®

## Automatic Handpiece Maintenance System

- Maximize the life span of handpieces
- Economical, efficient and time-saving
- Filtered air used to rotate the handpieces for safety
- Care for any brand of handpieces
- Disposable mist filtration system
- Detachable door for easy cleaning
- Easy replenishment of solution tank



Handpieces are one of the significant investments in the dental office. It is crucial to maintain their effectiveness, so the practitioner will have the cutting efficiency and durability that they need. CareMaster is a solution for the dynamic clinical needs.



### Individual controlled

**S:** Short Cycle for HS\* (Approx. 35s)  
**L:** Long Cycle for LS\* (Approx. 45s)  
**E:** Extra Long Cycle for SR\* (Approx. 50s)



### Air Regulator

Built-in air regulator to ensure proper maintenance.



### No mess

Spray mist is generally absorbed within the disposable filter sponges that can be changed.

\*HS: High Speed Handpiece, LS: Low Speed Handpiece, SR: Speed Reduction Handpiece

## Key Specifications

<b>Rated Voltage</b>	AC100V-240V, 50/60Hz
<b>Input Power</b>	35W
<b>Working Pressure</b>	0.35-0.40Mpa (3.5-4.0KGf/cm2)(50-60psi)
<b>Capacity</b>	350ml
<b>Dimensions</b>	282x300x370mm
<b>Kit Content</b>	CareMaster Unit, 1000ml CareMaster Solution Plus, Mist Filter Set (pack of 12 pieces), Oil Absorber Sheets (Pack of 20), Oil Replenish Funnel with Filter, Power Cord, Air Tube

## Order Information



2 Turbine Handpieces  
1 E-Type Attachment  
**CODE: HPL2010**



1 Turbine handpiece  
2 E-Type Attachments  
**CODE: HPL2011**



BeyesLube®  
Aerosole Lubricant  
**CODE: LC2010**



**Beyes Dental Canada Inc.**  
23-595 Middlefield Rd  
Scarborough, ON M1V 3S2  
Canada  
**Tel: 1-855-603-1888**  
**Email: info@beyes.ca**  
**Website: www.beyes.ca**



**Lotus NL B.V.**  
Koningin Julianaplein 10, le Verd,  
2595AA, The Hague,  
Netherlands

### Your Dealer for Beyes Products



Scan the QR code  
for more information

Printed in Canada  
File #: ENB022R3  
Rev.3/09.16.24

© Beyes Canada Inc 2024. All Rights Reserved.  
Specifications are subject to change without notice.  
All brands are holders of their respective trademarks.  
All Beyes products are manufactured in strict compliance  
with Health Canada regulations ensuring quality.  
Manufacturing facility is ISO 13485 certified under the  
Medical Device Single Audit Program (MDSAP).