

# Accurator<sup>®</sup>

## Apex Locator

The new and advanced Beyes Accurator has been developed to accurately measure the length of canals in all conditions, including dry, wet and where blood is presented, with industry leading accuracy.

### Accurate Measurement

Instant accurate measurement of canal length.

### Autoclavable Accessories

Autoclavable file holders and lip clips to prevent contamination between patients.

### No Zero-Adjustment


It is not necessary to set the device on zero before measuring each individual canal.

### Automatic Calibration

No manual calibration necessary. The Automatic calibration ensures accuracy and eliminates the effect of changes in temperature, moisture, etc. Inside the canal, even during treatment.

### Large Color LCD

The large color LCD screen is easy to read, and provides a clear progressive display with high contrast.



## Accurator Series

CODE AL2010  
MODEL Accurator A1

CODE ALP002  
MODEL File holders

CODE ALP001  
MODEL Probe cord

CODE ALP003  
MODEL Lip clip(4pc)