

Non-optic

# Maxso® M200E Smart

It is the most economical air driven handpiece in the market with high quality materials and craftsmanship. It's designed to last, perform, and meet all the daily operative needs. It has been a great product for the dealers to compete with direct sale companies and capture the price-conscious customer.




## Single Spray

The precise spray to ensure the best cooling effect to the treatment site.



## Push Button

The simple operation of the push-button chucking system provides quick bur change.



## Stainless Steel

Advanced technologies and materials provide the longest lifetime.



## Noise Level

Effective frequency optimization ensures a pleasant and low turbine noise

## Technical Data

	Mini Head	Standard Head
Head Size (Width x Height)	10.8mm x 13.0mm	12.5mm x 13.1mm
Drive Air Pressure (PSI)	40-60 PSI	40-60 PSI
Power (WATT)	19	21
Speed (RPM)	400,000	350,000
Rotary Instruments	FG burs according to EN ISO 1797-1:1995	
Ø / max. length (MM)	1.6 / 25	1.6 / 25
Chuck Type	Push Button	Push Button
Water Spray	Single	Single

## Order Info.

Midwest 4 Hole Connection



**Mini Head**

Code **HP2105**

Model **M200E-M/M4**

**Standard Head**

Code **HP2106**

Model **M200E-S/M4**

Borden 2 Hole Connection



**Mini Head**

Code **HP2108**

Model **M200E-M/B2**

**Standard Head**

Code **HP2109**

Model **M200E-S/B2**



**Beyes**<sup>®</sup>

*The Next Standard*<sup>®</sup>

**Beyes Dental Canada Inc.**

30 Casebridge Court  
Toronto, Ontario, M1B 3M5  
Canada

Tel: 1-416-281-1888

Email: [info@beyes.ca](mailto:info@beyes.ca)

Website: [www.beyes.ca](http://www.beyes.ca)

Your Dealer for Beyes Products



Printed in Canada.  
B110EN.REBOW1719

© Beyes Canada Inc 2017. All Rights Reserved.  
Specifications are subject to change without notice.  
All brands are holders of their respective trademarks.