

## Maxso® M200 Smart

The Maxso series high speed handpiece is designed to meet expectations in modern dental procedures. High quality components and comfortable grip with Micro-tex finish provide superior handling that you desire. Furthermore, the well aimed triple water spray will provide excellent coolant effects. Maxso series comes with variety of connections; as a result, you simply attach the Maxso to your existing system.



### Triple Water Spray

Effectively cools the entire surface of the bur that avoids overheating and provides better patient comfort.





### Micro-Tex Coating

The latest surface coating technology creates a micro-textured and scratch proof to enhance handling and durability.



### Compatibility




Compatible with a great variety of coupling systems and adapts great to every existing treatment units.



## Technical Data

## Standard Head

## Mini Head

Head Size(Width x Height)	12.5mm x 13.1mm	10.8 x 13.0mm
Drive Air Pressure(PSI)	40-45 PSI	40-45 PSI
Power(WATT)	21	19
Speed(RPM)	350,000	400,000
Rotary Instruments	FG burs according to EN ISO 1797-1:1995	FG burs according to EN ISO 1797-1:1995
Ø / max. length(MM)	1.6 / 25	1.6 / 25
Chuck Type	Push Button	Push Button
Water Spray	Triple	Triple
Midwest 4 Hole Connection		
	CODE <b>HP2008</b> MODEL <b>M200-S/M4</b>	CODE <b>HP2009</b> MODEL <b>M200-M/M4</b>
Borden 2 Hole Connection		
	CODE <b>HP2012</b> MODEL <b>M200-S/B2</b>	CODE <b>HP2013</b> MODEL <b>M200-M/B2</b>
Beyes® PD Connection		
	CODE <b>HP6022</b> MODEL <b>M200-S/PD</b>	CODE <b>HP6012</b> MODEL <b>M200-M/PD</b>
NSK® Connection		
	CODE <b>HP2036</b> MODEL <b>M200-S/N</b>	CODE <b>HP2037</b> MODEL <b>M200-M/N</b>

**Beyes®**  
The Next Standard

Beyes Dental Canada Inc.  
595 Middlefield Road  
Toronto, Ontario, M1V 3S2  
Canada  
Tel: 1-416-281-1888  
Email: [info@beyes.ca](mailto:info@beyes.ca)  
Website: [www.beyes.ca](http://www.beyes.ca)

## Your Dealer for Beyes Products

Scan the QR code for more information



Printed in Canada  
File #: ENB030  
Rev1/30.09.20

© Beyes Canada Inc 2020. All Rights Reserved.  
Specifications are subject to change without notice.  
All brands are holders of their respective trademarks.  
All Beyes products are manufactured in strict compliance with Health Canada regulations ensuring quality.  
Manufacturing facility is ISO 13485 certified under Canadian Medical Devices Conformity Assessment System (CMDACS).