


# AirLight<sup>®</sup>

## M800 Surgical 45

The AirLight M800 Surgical 45 provides excellent accessibility and visibility to make stress-free treatment a reality. It features 22W of cutting power for faster preparation and a rear exhaust to help increase patient safety. With its uniquely angled head, enhanced access is simple. Experience for yourself the difference between ordinary and extraordinary.




### Accessibility

45° Head design for perfect access to the molars.



### Visibility

Stream water (not mist) for perfect bur cooling and visibility.



### Hygiene

Rear air exhaust to prevent contamination into cutting area.




### Power

Strong chuck holding force to prevent slippage and improve power with 22W torque.

Generator+LED

## Technical Data

Colour temperature(K)	5,500
Brightness(LUX)	25,500
Head Size(Width x Height)	10.8mm x 13.0mm
Drive Air Pressure(PSI)	40-60 PSI
Power(WATT)	22
Speed(RPM)	410,000
Rotary Instruments	FG burs according to EN ISO 1797-1:1995
Ø / max. length(MM)	1.6 / 25
Chuck Type	Push Button
Water Port	Triple Water Jet



## Order Info.

Midwest 4 Hole Connection



Code **HP3048**  
Model **M800-45/M4**

Borden 2 Hole Connection



Code **HP3052**  
Model **M800-45/B2**

Beyes® PD Connection



Code **HP6038**  
Model **M800-45/PD**

Kavo® Connection



Code **HP3037**  
Model **M800-45/K**

NSK® Connection



Code **HP3044**  
Model **M800-45/N**

GLEED

**Beyes**  
CANADA

Beyes Dental Canada Inc.  
23-595 Middlefield Rd  
Scarborough, ON M1V 3S2  
Canada  
Tel: 1-855-603-1888  
Email: info@beyes.ca  
Website: www.beyes.ca



**Lotus NL B.V.**  
Koningin Julianaplein 10, Ie Verd,  
2595AA, The Hague,  
Netherlands

Your Dealer for Beyes Products



Scan the QR code  
for more information

Printed in Canada  
File #: ENB054R1  
Rev.1/09.16.24

© Beyes Canada Inc 2024. All Rights Reserved.  
Specifications are subject to change without notice.  
All brands are holders of their respective trademarks.  
All Beyes products are manufactured in strict compliance  
with Health Canada regulations ensuring quality.  
Manufacturing facility is ISO 13485 certified under the  
Medical Device Single Audit Program (MDSAP).