

Non-optic

# Maxso<sup>®</sup>

## M200 Surgical 45

It provides excellent accessibility, and visibility to make stress-free treatment a reality.

Surgi45 features 20W of cutting power for faster preparation and rear exhaust to help increase patient safety. With its uniquely angled head, enhanced access is simple. Experience for yourself, the difference between ordinary and extraordinary.



### Accessibility

45° Head design for perfect access to the molars.



### Visibility

Stream water (not mist) for perfect bur cooling and visibility.



### Hygiene

Rear air exhaust to prevent contamination into cutting area.



### Power

Strong chuck holding force to prevent slippage and improve power with 20W torque.

Non-optic

## Key Specifications

Head Size(Width x Height)	10.8mm x 13.0mm
Head Angle	45°
Drive Air Pressure(PSI)	40-60 PSI
Power(WATT)	20
Speed(RPM)	400,000
Ø / max. length(MM)	1.6 / 25
Chuck Type	Push Button
Water Spray	Single
Rotary Instruments	FG burs according to EN ISO 1797-1

## Order Information

Midwest 4 Hole Connection



Code **HP2010**  
Model **M200-45/M4**

Borden 2 Hole Connection



Code **HP2014**  
Model **M200-45/B2**

Beyes® QD Connection



Code **HP2018**  
Model **M200-45/QD**

**Beyes**<sup>®</sup>  
The Next Standard

Beyes Dental Canada Inc.  
595 Middlefield Road  
Toronto, Ontario, M1V 3S2  
Canada  
Tel: 1-416-281-1888  
Email: info@beyes.ca  
Website: www.beyes.ca

Your Dealer for Beyes Products

Scan the QR code for more information on M200 Surgical 45



Printed in Canada  
File #: ENB008  
Rev.2/30.09.20

© Beyes Canada Inc. 2020. All Rights Reserved.  
Specifications are subject to change without notice.  
All brands are holders of their respective trademarks.  
All Beyes products are manufactured in strict compliance with Health Canada regulations ensuring quality.  
Manufacturing facility is ISO 13485 certified under Canadian Medical Devices Conformity Assessment System (CMDCAS).