

Maxso[®] EP12-NS

Hygiene Handpieces

Beyes EP12-NS hygiene handpiece is a 4:1 reduction straight nose cone designed to be used with disposable prophy angle for polishing. The EP12-NS offers the highest quality ceramic bearings for long-lasting smooth operation and a slim, lightweight, and compact handpiece design for ultimate user comfort.

- Beyes EP12-NS will fit all standard ISO E-Type air motor or electric micromotor
- 4:1 reduction straight nose
- Ceramic bearings
- Amazingly quiet operation
- Autoclavable
- 1 Year warranty

Order Info

Code: SL1002
Model: EP12-NS



For Disposable
Prophy Angles



Key Specification

Transmission ratio	4:1
Connection	ISO E-Type
For use with	Disposable prophy angles
Speed type	Low speed
Maximum speed	10,000 rpm
Optic	Non optic
Overall length	64 mm
Shape	Straight
Weight	57 g
Sterilization	Autoclavable up to 135°C
Warranty	1 year warranty

Beyes
CANADA

Beyes Dental Canada Inc.
23-595 Middlefield Rd
Scarborough, ON M1V 3S2
Canada
Tel: 1-855-603-1888
Email: info@beyes.ca
Website: www.beyes.ca



Lotus NL B.V.
Koningin Julianaplein 10, le Verd,
2595AA, The Hague,
Netherlands

Your Dealer for Beyes Products



Scan the QR code
for more information

Printed in Canada
File #: ENB012R2
Rev.2/09.16.24

© Beyes Canada Inc 2024. All Rights Reserved.
Specifications are subject to change without notice.
All brands are holders of their respective trademarks.
All Beyes products are manufactured in strict compliance
with Health Canada regulations ensuring quality.
Manufacturing facility is ISO 13485 certified under the
Medical Device Single Audit Program (MDSAP).