

Maxso®
X200

X

traordinary



POWERFUL QUIET SAFE INNOVATIVE

Key Features

27W Torque

Increases efficiency and minimizes stall for the best precision cutting.

X Ball Bearings

Newly engineered ball bearings with special caging material that lasts four times longer than other ball bearings.

Quattro Spray

Four specially-positioned water ports guaranteeing optimal cooling.

Fiber Optic+

Provides a brighter and larger light pattern enabling clinicians with a clearer operative field.

IS Technology

Provides instant stop capability which eliminates backflow and noise when ceasing operation.

Micro-Tex

A slip-proof coating improves hand grip while ensuring safety and comfort.

The well-designed and precisely-manufactured Maxso X200X delivers an extraordinary torque of **27 watts** ensuring accurate cutting. This handpiece is equipped with specially-designed X ball bearings for quieter sound and longer life. It also possesses instant-stop capability (**IS Tech**) effectively reducing noise and eliminating backflow. In addition to four specially-positioned water ports guaranteeing optimal cooling, it has an all-new **Fiber Optic+** that provides a larger and brighter light pattern, enabling clinicians with a clearer operative field. The Maxso X200X also comes with an extended **25 month warranty**.



Beyes® Handpieces

Model	Maxso® X200X	AirLight® M800	AirLight® M800P	Maxso® X200S	Maxso® M200S	Maxso® M200E
RPM	430,000	400,000	400,000	400,000	400,000	400,000
Power	27	23	23	23	21	21
Spray Ports	4	3	3	3	3	1
Illumination	Fiber Optic+	LED	LED	Fiber Optic	—	—
Instant Stop	✓	—	—	—	—	—
Material	Stainless Steel	Special Alloy	Titanium	Special Alloy	Special Alloy	Special Alloy
Warranty	25	12	12	12	12	6



Key Specifications

Power (Watts)	27
Drive Air Pressure (PSI)	40-50
Speed (RPM)	430,000 (10%±)
Water Spray	Quattro
Colour Temperature (K)	5,500
Brightness (LUX)	25,500
Ø / Max. Length (MM)	1.6 / 25
Chuck Type	Push Button
Head Size - W x H (MM)	10.8 x 13.0
Rotary Instruments	FG burs according to EN ISO 1797-1

Order Information

Kavo® Backend



Code: **HP2060X**
Model: **X200X-M/K**

NSK® Backend



Code: **HP2066X**
Model: **X200X-M/N**

Beyes®

Partnership for Success

Beyes Dental Canada Inc.
23-595 Middlefield Rd
Scarborough, ON M1V 3S2
Canada

Tel: 1-855-603-1888

Email: info@beyes.ca

Website: www.beyes.ca



Lotus NL B.V.
Koningin Julianaplein 10, 1e Verd,
2595AA, The Hague,
Netherlands

Your Dealer for Beyes Products



Scan the QR code for
more products

Printed in Canada
File #: ENB010
Rev.3 /09.17.24

© Beyes Canada Inc 2024. All Rights Reserved.
Specifications are subject to change without notice.
All brands are holders of their respective trademarks.
All Beyes products are manufactured in strict compliance
with Health Canada regulations ensuring quality.
Manufacturing facility is ISO 13485 certified under the
Medical Device Single Audit Program (MDSAP).