

Maxso®
X200

X

traordinary



POWERFUL QUIET SAFE INNOVATIVE

Key Features

X Ball Bearings

Newly engineered ball bearings with special caging material that lasts four times longer than other ball bearings.

Quattro Spray

Four specially-positioned water ports guaranteeing optimal cooling.

27W Torque

Increases efficiency and minimizes stall for the best precision cutting.

IS Technology

Provides instant stop capability which eliminates backflow and noise when ceasing operation.

Fiber Optic+

Provides a brighter and larger light pattern enabling clinicians with a clearer operative field.

Micro-Tex

A slip-proof coating improves hand grip while ensuring safety and comfort.

The well-designed and precisely-manufactured Maxso X200X delivers an extraordinary torque of **27 watts** ensuring accurate cutting. This handpiece is equipped with specially-designed X ball bearings for quieter sound and longer life. It also possesses instant-stop capability (**IS Tech**) effectively reducing noise and eliminating backflow. In addition to four specially-positioned water ports guaranteeing optimal cooling, it has an all-new **Fiber Optic+** that provides a larger and brighter light pattern, enabling clinicians with a clearer operative field. The Maxso X200X also comes with an extended **25 month warranty**.



X200X

Beyes® Handpieces

Model	Maxso® X200X	AirLight® M800	AirLight® M800P	Maxso® X200S	Maxso® M200S	Maxso® M200E
RPM	430,000	400,000	400,000	400,000	400,000	400,000
Power	27	23	23	23	21	21
Spray Ports	4	3	3	3	3	1
Illumination	Fiber Optic+	LED	LED	Fiber Optic	—	—
Instant Stop	✓	—	—	—	—	—
Material	Stainless Steel	Special Alloy	Titanium	Special Alloy	Special Alloy	Special Alloy
Warranty	25	12	12	12	12	6



Key Specifications

Power (Watts)	27
Drive Air Pressure (PSI)	40-50
Speed (RPM)	430,000 (10%±)
Water Spray	Quattro
Colour Temperature (K)	5,500
Brightness (LUX)	25,500
Ø / Max. Length (MM)	1.6 / 25
Chuck Type	Push Button
Head Size - W x H (MM)	10.8 x 13.0
Rotary Instruments	FG burs according to EN ISO 1797-1

Order Information

Kavo® Backend



Code: **HP2060X**
Model: **X200X-M/K**

NSK® Backend



Code: **HP2066X**
Model: **X200X-M/N**

Beyes®
Partnership for Success

Beyes Dental Canada Inc.
23-595 Middlefield Rd.
Toronto, Ontario, M1V 3S2
Canada
Tel: 1-416-281-1888
Email: info@beyes.ca
Website: www.beyes.ca



Lotus NL B.V.
Koningin Julianaplein 10, 1e Verd,
2595AA, The Hague,
Netherlands
Tel: +31645171879 (English)
+31626669008 (Dutch)
Website: www.beyes.ca

Your Dealer for Beyes Products

Scan the QR code for
more products



Printed in Canada
File #: ENB010
Rev.2 /05.02.20

© Beyes Canada Inc 2018. All Rights Reserved.
Specifications are subject to change without notice.
All brands are holders of their respective trademarks.
All Beyes products are manufactured in strict compliance
with Health Canada regulations ensuring quality.
Manufacturing facility is ISO 13485 certified under Canadian
Medical Devices Conformity Assessment System (CMDCAS).