

**Summary of ADA professional Product Review**  
**ADA is abbreviation of American Dental Association**  
**Key Features ONLY**  
**Volume 10. ISSUE 1. 2015**

Page #	Specifications	W&H/A-Dec	Beyes Canada	NSK/Brasseler	KaVo	Lares	Dentsply	Dentsply	Star Dental
		Synea TA-97 LED+	AirLight Mini	Mini KV/Pico	M8900L	557 Ultralite	Tradition	Stylus Lite	Solara QT-L
5	Price (\$US Dollar)	\$1,650	\$698	\$845	\$1,550	\$770	N/A	N/A	N/A
8	Maximum Diameter of Head (mm)	10	10.9	8.7	12.5	10.2	11.2	10.6	10.6
8	Overall Length (mm)	120	120	128	127	117	126	125	118
9	Speed (Rotation Per Min/RPM)	409,465	419,376	387,915	355,031	492,776	370,607	412,042	397,950
14	Sound Level (Decibels)	73	72.3	69.7	69.2	73.7	70.6	72.4	73.4
20	Power/Torque (Watt)	14.5	17.1	7	21.8	13.4	11.3	11.8	14.4
22	Illumination	LED	Direct-LED	Fiber Optic	Fiber Optic	Fiber Optic	Fiber Optic	Fiber Optic	Fiber Optic

Original Source: [http://www.ada.org/~media/ADA/Publications/PPR/2015%20PPR%20Issues/PPR\\_VOL10\\_ISS1.ashx](http://www.ada.org/~media/ADA/Publications/PPR/2015%20PPR%20Issues/PPR_VOL10_ISS1.ashx)