

Maxso® E800

Electric Handpiece System



Refined Trimming
Faster Tooth Preparation
More Powerful Crown-splitting

Beyes® has made upgrading to the newest electric technology easy with the Maxso® E800. Equipped with a touch screen interface and 3 programmable preset speeds, this brushless electric micromotor system is designed with ultimate user-friendliness in mind. With customizable settings, the system is the perfect solution for all high-speed and low-speed applications. Its quiet motor delivers uncompromised power – featuring advanced rotational control technology, resulting in only the smoothest and steadiest operation. Integrate the Maxso® E800 into your practice with ease and start experiencing the advantages of electric today.

The E800 Advantage:

- **Powerful:** 5.0 Ncm of torque
- **Steady Performance:** Speed of rotations is controlled between 2,000 and 40,000min-1
- **Smooth Operation:** Torque ripple and heat reduction technology eliminates vibrations
- **Convenient & Durable:** Smooth and easy to clean surface with autoclavable brushless motor
- **3 Programmable Memory Settings:** Easily switch between speed presets for more efficient workflow
- **New Touch Screen Interface:** Easily and quickly navigate through menu settings on a color LCD display panel

Maxso® E800P



Key Specification

Motor speed range	2000-40000rpm
Maximum torque	5.0 N.cm
Bur offset	<1.3um
Compatible contra-angle	16:1, 1:1, 1:5(Standard)
Motor size	Φ22x76.7mm
Motor weight	70g
Screen Dimensions	W70 x H45 (mm)
Noise Level	<60dB
Warranty	25 months

Order Information

Portable Electric Handpiece System

Code: MT8001
Model: Maxso E800P

Built-in Electric Handpiece System

Code: MT8002
Model: Maxso E800B

MC Motor Only, Brushless

Code: MT800MC
Model: Maxso E800 MC



Beyes Dental Canada Inc.
23-595 Middlefield Road
Toronto, Ontario, M1B 3M5
Canada
Tel: 1-855-603-1888
Email: info@beyes.ca
Web: www.beyes.ca

Lotus NL B.V.
Koningin Julianaplein 10, le Verd,
2595AA, The Hague,
Netherlands
Tel: +31645171879 (English)
+31626669008 (Dutch)
Email: info@beyes.ca

Your Dealer for Beyes Products

Scan the QR code
for more products



Printed in Canada
File #: ENB041
Rev.1/07.07.21

© Beyes Canada Inc 2021. All Rights Reserved.
Specifications are subject to change without notice.
All brands are holders of their respective trademarks.
All Beyes products are manufactured in strict compliance
with Health Canada regulations ensuring quality.
Manufacturing facility is ISO 13485 certified under Canadian
Medical Devices Conformity Assessment System (CMDCAS).