

Comfortek®

Dental Cart System

Comfortek Dental Cart System from Beyes is designed ergonomically with ultimate user comfort in mind. This complete mobile dental cart is height adjustable and has foot pedal control of water switch and suction to make your work quicker and easier. Along with a smooth tabletop design for easy cleaning, there is plenty of space for all necessary instruments.

The complete system includes Beyes AirLight M800 Plus high speed handpiece - a direct LED handpiece that produces brighter illumination than fiber-optics, without needing fiber-optic tubing. However, we understand that there may be varying needs across different dental facilities. As such, Comfortek Dental Cart System is sold as a fully customizable unit.

Comfortek Dental Cart System also features an eco-friendly oil-free compressor that decreases the risk of contamination and requires less maintenance from your staff.



Features

- Suction with Foot Pedal Control
- 110V Oil-free Compressor
- Air/Water Syringe with 360° Tip Rotation
- Autoclavable LED Piezo Scaler
- Height Adjustable from 33.5" to 44.5"
- Top Shelf with Rounded Edges
- 2 Water Bottles with Toggle Switch
- Safety Toggle Switch to Depressurize Water System
- Foot Control (Water Control, Suction Control, Unit Control)
- Separate On/Off Switch For Compressor
- Water and Air Flow Adjustments for Individual Handpieces
- 2.5" Castors with Brakes



Key Specifications

Top Shelf Measurement	15.75"(L) x 13"(W)
Shelf Base Measurement	17.75"(L) x 17.75"(W) x 6.5"(H)
Base Measurement	17.75"(L) x 17.75"(W)
Adjustable Height	33.94"(min) - 45.28" (max)
Weight	56kg
Capacity	15kg
AMPS	6.0A
Power Source	115VAC
Operating Frequency	50~60HZ
Tank Capacity	5L(1.32gal)
Air Pressure	7.5bar(108psi)
Warranty	3 years

Order Information

Code: MB1006

Model: C1

Description: Comfortek Dental Cart System, With Air Compressor, Suction

Code: MB1007

Model: C2

Description: Comfortek Dental Cart System, With Air Compressor

Code: MB1008

Model: C3

Description: Comfortek Dental Cart System, Without Air Compressor, With Suction

Code: MB1009

Model: C4

Description: Comfortek Dental Cart System, Without Air Compressor

Components

(For complete system only)

- ▶ **Air/Water Syringe**
- ▶ **LED Piezo Scaler**
 - 5 Scaling/Perio Tips Included
 - Torque Wrench
- ▶ **LED High-speed Handpiece**
 - AirLight M800 Plus Handpiece
 - Beyes PD Coupler
- ▶ **Low Speed Handpiece**
 - 20K Air Motor
 - 4:1 Straight Nose Cone for DPA
 - Set of 200 Disposable Prophy Angles
- ▶ **Suction**
 - Saliva Injector Tip
 - Suction Curved Tip



Beyes
CANADA

Beyes Dental Canada Inc.
23-595 Middlefield Rd
Scarborough, ON M1V 3S2
Canada
Tel: 1-855-603-1888
Email: info@beyes.ca
Website: www.beyes.ca



Lotus NL B.V.
Koningin Julianaplein 10, le Verd,
2595AA, The Hague,
Netherlands

Your Dealer for Beyes Products



Scan the QR code
for more information

Printed in Canada
File #: ENB006R4
Rev.3/09.26.24

© Beyes Canada Inc 2024. All Rights Reserved.
Specifications are subject to change without notice.
All brands are holders of their respective trademarks.
All Beyes products are manufactured in strict compliance
with Health Canada regulations ensuring quality.
Manufacturing facility is ISO 13485 certified under the
Medical Device Single Audit Program (MDSAP).