

Sensibles™

Universal Sensor Holder

Make the Sensible choice!

Engineered for adaptability. Designed for simplicity.

- All-in-one solution for the perfect image
- Simple and secure adjustment
- Perfect for Anterior, Posterior, and Bite-Wing images
- Autoclavable and reusable
- Durable and 100% latex-free



Starter Kit Includes

- 2 x Bite Block Large Size
- 2 x Bite Block Medium Size
- 1 x Posterior Aiming Ring
- 1 x Anterior/Bite Wing Aiming Ring

- 1 x Positioning Arm Posterior
- 1 x Positioning Arm Bite Wing
- 1 x Positioning Arm Anterior



Horizontal Imaging

The Medium Sensible holds most size 2, 1 and 0 sensors horizontally.



Vertical Imaging

The Large (two-tone) Sensible holds most size 2 and 1 sensors vertically.



Bite-Wing Imaging

For Bite Wing imaging, simply slide Bite Block to mid-line of sensor.


USING SENSIBLES IS AS EASY AS 1-2-3:



1 Attach Positioning Bar



2 Add Aiming Ring



3 Snap In Sensor



For Bite-Wings, adjust vertical position

Order Information

Order Code	Model / Description	Order Code	Model / Description
DS1103	Universal Starter Kit	DS1107	Bite Block Large Size Refill (12/pk)
DS1104	Positioning Arm Posterior (1/pk)	DS1108	Bite Block Medium Size Refill (12/pk)
DS1105	Positioning Arm Bite Wing (1/pk)	DS1109	Posterior Aiming Ring (1/pk)
DS1106	Positioning Arm Anterior (1/pk)	DS1110	Anterior/Bite Wing Aiming Ring (1/pk)

