

BeyesLube®

Advanced Auto-Lube



FOR USE WITH AUTOMATED HANDPIECE LUBRICATION STATIONS

BeyesLube® Advanced Auto-Lube provides superior autoclave-safe lubrication and eliminates costly repairs. It extends turbine life by reducing friction and wear, restoring bearing cages and o-rings, and preventing buildup that causes sluggish handpieces. For use in all automatic dental handpiece lubrication and cleaning stations including CareMaster (Beyes), Handpiece Spa (Tuttnauer), Assistina (A-dec/W&H), and Care 3 (NSK). BeyesLube® is rated H-1 Food Grade, non-flammable, non-hazardous, and meets and/or exceeds all published manufacturers' specifications.

1000ml
(33.8 fl oz)



Order Code: HPLP012

Beyes[®]
The Next Standard[®]



WARNING

Hazardous Statement:

General:

If medical advice is needed, have product label at hand. Keep out of reach of children. Read label before use.

Precautionary Statements:

Keep away from hot surfaces, sparks, open flames and other ignition sources. No smoking. Do not pierce or burn, even after use. Protect from sunlight. Do not expose to temperatures exceeding 50°C / 122°F.

First Aid:

IF ON SKIN: Wash with plenty of water. If skin irritation occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse.

IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

Ingestion Information:

Immediately call a POISON CENTER/doctor. Do not induce vomiting. If vomiting occurs naturally, lie on your side, in the recovery position. Do not give anything.



Beyes Dental Canada Inc.
23-595 Middlefield Rd
Scarborough, ON M1V3S2
Canada
Tel: 1-416-281-1888
Email: info@beyes.ca
Web: www.beyes.ca



Lotus NL B.V.
Koningin Julianaplein 10, le Verd,
2595AA, The Hague,
Netherlands

Lot #

