

AirLight® M800 Plus

LED High Speed Handpiece

- 23W of Torque
- 410,000 RPM
- Ceramic Ball Bearings
- Handles up to 60 PSI of air pressure
- Titanium Body
- 1 Year Warranty

Backend	Code
4 Hole	HP3008P
Beyes®	HP6018P
W&H®	HP3022P
Star®	HP3054P
NSK®	HP3043P
KAVO®	HP3036P

Reg. \$548/ea
Net: \$365³³/ea

BUY 2, GET 1 FREE!



Direct LED Without Fiber Optics!
Just Brilliantly Bright Illumination!

Built-in
Generator

BUY 2, GET 1 FREE!

PD Coupler
Non Optic, 4 Holes

Code: HP6008

Reg. \$152/ea **Net: \$101.33/ea**

Maxso® E800

Electric Micromotor System Package

- Powerful: 5.0 Ncm of torque
- Steady Performance:
Speed of rotations is controlled between 2,000 and 40,000min-1
- Smooth Operation:
Torque ripple and heat reduction technology eliminates vibrations
- Convenient & Durable:
Smooth and easy to clean surface with autoclavable brushless motor
- 3 Programmable Memory Settings:
Easily switch between speed presets for more efficient workflow
- New Touch Screen Interface:
Easily and quickly navigate through menu settings on a color LCD display panel
- Warranty: 25 months

MT8001KIT1:
E800P+One X99L Attachment
Reg. \$2,398 **Sale: \$2,075**

MT8001KIT2:
E800P+Two X99L Attachment
Reg. \$3,118 **Sale: \$2,699**

Maxso® X99L

1:5 Contra Angle

- Anti-Heat System
- Ceramic Ball Bearing

- ISO E-Type Connection
- Quattro Spray
- Fiber-optic
- 1 Year Warranty

BUY 2, GET 1 FREE!

Code: SL3025
Reg. \$992/ea
Net: \$661.33/ea

Maxso® Low Speed

E-Type Connection

BUY 3, GET 1 FREE!



1 Year Warranty **Mix and Match**

- A. SL2026** 1:1 Push Button (CA Burs)
B. SL2028 1:1 Latch Type (CA Burs)
C. SL2034 1:1 Screw-On
D. MT2010 Air Motor, 4 Hole, 20,000RPM
E. ST2010 1:1 Straight Nose Cone
F. SL1002 4:1 Nose Cone, For DPA

Reg
\$176/ea
\$147/ea
\$139/ea
\$298/ea
\$114/ea
\$278/ea

Net
\$132/ea
\$110²⁵/ea
\$104²⁵/ea
\$223⁵⁰/ea
\$85⁵⁰/ea
\$208⁵⁰/ea

Mix & Match any 3. Promotional goods must be of equal or lesser value than the lowest priced item purchased.

AirLight® M800X

Generator + LED High Speed

- 27W of Torque
- 430,000 RPM
- Instant Stop Capability
- X-Ball Bearings
- Extra Quiet
- Micro-Tex Slip-Proof Coating
- 25 Months Warranty

Backend	Code
Beyes®	HP6018X
Star®	HP3054X
NSK®	HP3043X
KAVO®	HP3036X

Reg. \$1,275/ea
Net: \$637⁵⁰/ea

BUY 1, GET 1 FREE!



Built-in
Generator

Airlight® M800 Surgical 45

LED High Speed Handpiece

- 22W of Torque
- 45 Degree Angled Head
- 410,000 RPM
- Rear Air Exhaust
- Direct LED
- Stream Water
- 1 Year Warranty

Backend	Code
4 Hole	HP3048
Beyes®	HP6038
NSK®	HP3044
KAVO®	HP3037
Star®	HP3055

Reg. \$598/ea
Net: \$448⁵⁰/ea

BUY 3, GET 1 FREE!



Built-in
Generator

Canaview®

Intraoral Camera

- True Auto-Focus
- Anti-Fog Design
- Responsive Capture Button
- 5 MP High Resolution Images
- Intuitive Camera Lighting System
- No Manual Adjustments Required
- Direct USB Connection to Computer
- Super Light-Weight
- Complimentary Image Management Software

Code: EV2060

Reg. \$904/ea **Net: \$452/ea**

BUY 1, GET 1 FREE!



Slimax-C Plus®

LED Curing Light System

- Built-In Digital Radiometer
- LED Display Panel
- Cordless, Portable Design
- Cures 850 Times Per Charge
- Up to 2000 mW/cm ± 10%
- Rechargeable Samsung Lithium-Ion Battery
- Solid Aluminum Body
- 1 Year Warranty

Color	Code
Silver	SM1003P-S
Blue	SM1003P-B
Pink	SM1003P-Pi

Reg. \$660/ea
Net: \$330/ea

BUY 1, GET 1 FREE!



FREE GOODS REDEMPTION: 1. Go to "http://beyes.ca/en/redeems/new" 2. Fill out the form and submit (Promo code: **USQ124**)

**Free goods must be of equal or lesser value and cannot be combined with any other offers. All offers will require a proof of purchase dated from January 1st - March 31st, 2024 and must be reflected on one invoice. Offers must be redeemed no later than April 15, 2024. Please allow 4 to 8 weeks for promotion to be processed. All brands are holders of their respective trademarks.