

Comfortek® X1

Dental Chair with Delivery Unit

Includes:

- Stitched upholstery design
- Heavy-duty aluminum backrest
- Double post design
- Flexible delivery system with control panel
- 3 programmable memory positions
- One 6-pin fiber-optic tubing
- Built-in Piezo scaler
- Removable and rotatable stainless steel water outlet
- 4 instrument position on premium swiveling assistant side with premium control panel
- Main switch for all (air, water and electricity)
- 4-way foot switch
- Dentist and assistant stool

Now Upgraded with

- M7 LED Light with Composite Safe Mode
- Extended 8 years warranty

Reg. \$12,999

Sale: \$9,999



AirLight® M800 Plus

LED High Speed Handpiece

- 28W of Torque
- 410,000 RPM
- Ceramic Ball Bearings
- Handles up to 60 PSI of air pressure
- Titanium Body
- 12 Months Warranty

Backend	Code
4 Hole	(HP3008P)
Beye®	(HP6018P)
W&H®	(HP3022P)
Star®	(HP3054P)
NSK®	(HP3043P)
KAVO®	(HP3036P)

Reg. \$625/ea

Net: \$416⁶⁷/ea

Triple Spray

Direct-LED

Built-in Generator



BUY 2, GET 1 FREE

Airlight® M800 Surgical 45

LED High Speed Handpiece

- 25W of Torque
- 45 Degree Angled Head
- 410,000 RPM
- Rear Air Exhaust
- Direct LED
- Stream Water
- 12 Months Warranty

Backend	Code
4 Hole	(HP3048)
Beyes®	(HP6038)
NSK®	(HP3044)
KAVO®	(HP3037)
Star®	(HP3055)

Reg. \$680/ea

Net: \$453³³/ea

Triple Spray

Direct-LED

Built-in Generator



BUY 2, GET 1 FREE

AirLight® M800X

LED High Speed Handpiece

- 35W of Torque
- 430,000 RPM
- Extra Quiet
- 25 Months Warranty

Torque 35w

Backend	Code
Beyes®	(HP6018X)
Star®	(HP3054X)
NSK®	(HP3043X)
KAVO®	(HP3036X)

Reg. \$1,450/ea **Net: \$725/ea**

Limitless PSI Technology

Experience the freedom from operating pressure control

X Ball Bearings

Lasts 4x longer than other ball bearings

Instant Stop Technology

Eliminates backflow and prevents cross-contamination



BUY 1, GET 1 FREE

Slimax-C Plus®

LED Curing Light System

- Built-In Digital Radiometer
- LED Display Panel
- Cordless, Portable Design
- Cures 850 Times Per Charge
- Up to 2000 mW/cm ± 10%
- Rechargeable Samsung Lithium-Ion Battery
- Solid Aluminum Body
- 12 Months Warranty

Ultra-fast 5-Second Cure Time!

Color	Code
Silver	(SM1003P-S)
Blue	(SM1003P-B)
Pink	(SM1003P-PI)

Reg. \$625/ea **Net: \$312⁵⁰/ea**



BUY 1, GET 1 FREE

Maxso® X200X

Fiber-Optic High Speed Handpiece

- 35W of Torque
- 430,000 RPM
- Instant Stop Capability
- Quattro Spray
- Titanium Body
- 25 Month Warranty

Torque 35w

Backend	Code
NSK®	HP2066X
KAVO®	HP2060X

Reg. \$1,200/ea **Net: \$600/ea**



BUY 1, GET 1 FREE

Canaview®

Intraoral Camera

- 5MP High resolution images
- True auto-focus and anti-fog lens
- Plug-and-play with direct USB connection
- Intuitive camera lighting system
- Super lightweight
- Compatible with all major imaging software including Dentrax, Dexis, Eaglesoft, Curve Hero, Open Dental ...and more!
- 12 Months Warranty



Reg. \$999/ea **Net: \$499⁵⁰/ea** Code: (EV2060)

BUY 2, GET 2 FREE

FREE GOODS REDEMPTION: 1. Go to "http://beyesdental.com/en/redeems/new" **2.** Fill out the form and submit (Promo code: **USQ126**)

*Free goods must be of equal or lesser value and cannot be combined with any other offers. All offers will require a proof of purchase dated from January 1 - March 31, 2026 and must be reflected on one invoice. Offers must be redeemed no later than April 15, 2026. Please allow 4 to 8 weeks for promotion to be processed. All brands are holders of their respective trademarks. All promotions are subject to change or cancellation without notice. Beyes and its authorized distributors reserve the right to modify offers, correct errors and limit quantities at any time; in the event of a discrepancy, the price and terms confirmed at the point of sale will prevail.

Maxso® E800

Electric Handpiece System

- **Powerful:** 5.0 Ncm of torque
- **Steady Performance:** Speed of rotations is controlled between 2,000 and 40,000min-1
- **Smooth Operation:** Torque ripple and heat reduction technology eliminates vibrations
- **Convenient & Durable:** Smooth and easy to clean surface with autoclavable brushless motor
- **3 Programmable Memory Settings:** Easily switch between speed presets for more efficient workflow
- **New Touch Screen Interface:** Easily and quickly navigate through menu settings on a color LCD display panel
- **Warranty:** 25 Months

Code: MT8001
Reg. \$1,899



(1) E800P + (1) X99L
Code: MT8001KIT1
Reg. \$2,700
Sale: \$2,090



(1) E800P + (2) X99L
Code: MT8001KIT2
Reg. \$3,500
Sale: \$2,599



BUY 2, GET 1 FREE

Maxso® X99L

1:5 Contra Angle

- Anti-Heat System
- Quattro Spray
- Ceramic Ball Bearing
- Fiber-optic
- ISO E-Type Connection
- 12 Months Warranty



Code: SL3025

Reg. \$1,090/ea **Net: \$726^{.67}/ea**

Maxso® Electric Low speed

Internal Spray

BUY 2, GET 1 FREE

Mix & Match

12 Months Warranty



	Reg	Net
A. SL3068 M29L-MH01, Push Button (CA Burs), Fiber-Optic	\$550	\$366.67
SL3066 M29-MH01, Push Button (CA Burs), Non-Optic	\$300	\$200
B. ST2008 S20A-IS, 1:1 Straight Nose Cone, Non-Optic	\$305	\$203.33
ST2009 S20A-IS/E6, 1:1 Straight Nose Cone, Fiber-Optic	\$520	\$346.67

Mix & Match any 3. Promotional goods must be of equal or lesser value than the lowest priced item purchased.

EZ-Swivel®

Couplers

Mix & Match

BUY 2, GET 1 FREE

	Reg	Net
A. HP6008 Beyes® PD, 4 Holes, Non-Optic	\$170	\$113.33
B. HP3046 NSK® Type, 4 Holes, Non-Optic	\$285	\$190
HP2069 NSK® Type, 6 Holes, Fiber-Optic	\$285	\$190
C. HP2063 Kavo® Type, 4 Holes, Non-Optic	\$285	\$190
HP2064 Kavo® Type, 6 Holes, Fiber-Optic	\$285	\$190



Mix & Match any 2. Promotional goods must be of equal or lesser value than the lowest priced item purchased.

Duray® Intraoral Sensor Package

- Thinnest Sensor 4.4mm
- Durable and Reliable: Waterproof with IP68 protection and can be immersed in disinfectant liquids to prevent cross infections
- High Resolution Imaging: With a theoretical resolution of 25 lp/mm and through advanced CMOS technology, image quality is exceptional
- Compatible with Major Imaging Software: Apeteryx, Dexis, Dentrax Ascend, Eaglesoft, Open Dental, Suni, Schick, Tigerview, Vixwin, XDR, and more!
- 24 Months Warranty



Model	Code	Reg.	Sale
Size 0	(DS0001)	\$2,399	\$2,099
Size 1	(DS1001)	\$3,200	\$2,650
Size 1.5	(DS1501)	\$3,500	\$2,850
Size 2	(DS2001)	\$3,999	\$3,250
Size 1+2 Package	(DS3001)	\$6,999	\$5,600



Maxso® Low Speed

Non-spray

BUY 3, GET 1 FREE

Mix & Match



12 Months Warranty

	Reg	Net
A. SL2026 1:1 Push Button (CA Burs)	\$190	\$142.50
B. SL2028 1:1 Latch Type (CA Burs)	\$165	\$123.75
C. SL2034 1:1 Screw-On	\$150	\$112.50
D. ST2010 1:1 Straight Nose Cone	\$145	\$108.75
E. SL1002 4:1 Nose Cone, For DPA	\$310	\$232.50

Mix & Match any 3. Promotional goods must be of equal or lesser value than the lowest priced item purchased.

CareMaster®

Automatic Handpiece Maintenance System

Maintain Your Most Significant Investments in the Office!
Clean 3 Handpieces Per Cycle!

- Maximize the lifespan of handpieces
- Economical, efficient and time-saving
- Filtered air used to rotate the handpieces
- Detachable door for easy cleaning
- Provides care for all brands of handpieces
- 12 Months Warranty



* Requires Airline Connection
* Coupler Sold Separately

Code: HPL2010
2 Turbines & 1 E-Type Attachment
Reg. \$1,350/ea

**Buy Any Model of CareMaster
Get 1 AirLight® M800 Plus FREE**

Code: HPL2011
1 Turbines & 2 E-Type Attachments
Reg. \$1,475/ea




Duray® 2D, Ceph & 3D

Extraoral X-ray Systems

Advanced Imaging for Every Practice

- Quick Cephalometrics with Touchscreen Access
- 2D, Ceph, and Fully Upgradeable 3D Options
- Automatic Metal Artifact & Noise Reduction (AMAR & ANR)
- Fast & Easy Patient Positioning | Plug & Play Installation
- Smart Backlight | Image Mirroring | Advanced Sensor
- 5 Year Warranty Across All Models

Model	Code
Duray 2D Pan	(XR3020)
Duray 2D Ceph	(XR3021)
Duray 3D	(XR3030)
Duray 3D Ceph	(XR3031)



Contact us for
complete specs
and pricing

www.beyesdental.com
info@beyes.ca
1-855-603-1888

FREE GOODS REDEMPTION: 1.Go to "http://beyesdental.com/en/redeems/new" 2.Fill out the form and submit (Promo code: USQ126)

*Free goods must be of equal or lesser value and cannot be combined with any other offers. All offers will require a proof of purchase dated from January 1 - March 31, 2026 and must be reflected on one invoice. Offers must be redeemed no later than April 15, 2026. Please allow 4 to 8 weeks for promotion to be processed. All brands are holders of their respective trademarks. All promotions are subject to change or cancellation without notice. Beyes and its authorized distributors reserve the right to modify offers, correct errors and limit quantities at any time; in the event of a discrepancy, the price and terms confirmed at the point of sale will prevail.