

**AirLight® M800 Plus**

LED High Speed Handpiece

- 23W of Torque
- 410,000 RPM
- Ceramic Ball Bearings
- Handles up to 60 PSI of air pressure
- Titanium Body
- 1 Year Warranty

Backend	Code
4 Hole	HP3008P
Beyes®	HP6018P
W&H®	HP3022P
Star®	HP3054P
NSK®	HP3043P
KAVO®	HP3036P

Reg. \$522/ea  
**Net: \$348/ea**

**BUY 2, GET 1 FREE!**



**Direct LED Without Fiber Optics!  
Just Brilliantly Bright Illumination!**

**Built-in  
Generator**



**AirLight® M800X**

Generator + LED High Speed

- 27W of Torque
- 430,000 RPM
- Instant Stop Capability
- X-Ball Bearings
- Extra Quiet
- Micro-TEX Slip-Proof Coating
- 25 Months Warranty

Backend	Code
Beyes®	HP6018X
Star®	HP3054X
NSK®	HP3043X
KAVO®	HP3036X

Reg. \$1,250/ea **Net: \$625/ea**

**BUY 1, GET 1 FREE!**



**Built-in  
Generator**



**BUY 2, GET 1 FREE!**

**PD Coupler**  
Non Optic, 4 Holes



Code: HP6008

Reg. \$144/ea **Net: \$96/ea**

**Duray® Intraoral Sensor Package**

**HOT DEAL**

- Thinnest Sensor 4.4mm
- Durable and Reliable  
Waterproof with IP68 protection and can be immersed in disinfectant liquids to prevent cross infections
- High Resolution Imaging  
With a theoretical resolution of 25 lp/mm and through advanced CMOS technology, image quality is exceptional
- Easy To Use  
Simple and easy-to-operate software with one-click smart installation and offers plug-and-play simplicity with high-speed USB 2.0 connection
- 24 Months Warranty



**(DS3001) Duray Sensor Package**

Reg. \$5999 **Sale: \$4999**

**(DS1001) Size 1**

Reg. \$2800 **Sale: \$2350**

**(DS2001) Size 2**

Reg. \$3500 **Sale: \$2999**



**Maxso® E800**

Electric Micromotor System Package

- Powerful: 5.0 Ncm of torque
- Steady Performance:  
Speed of rotations is controlled between 2,000 and 40,000min-1
- Smooth Operation:  
Torque ripple and heat reduction technology eliminates vibrations
- Convenient & Durable:  
Smooth and easy to clean surface with autoclavable brushless motor
- 3 Programmable Memory Settings:  
Easily switch between speed presets for more efficient workflow
- New Touch Screen Interface:  
Easily and quickly navigate through menu settings on a color LCD display panel
- Warranty: 25 months



**MT8001KIT1:**  
E800P+One X99L Attachment  
Reg. \$2,300 **Sale: \$1,999**



**MT8001KIT2:**  
E800P+Two X99L Attachment  
Reg. \$2,998 **Sale: \$2,599**

**Maxso® X99L**

- 1:5 Contra Angle
- Anti-Heat System
- Ceramic Ball Bearing

- ISO E-Type Connection
- Quattro Spray
- Fiber-optic
- 1 Year Warranty



**BUY 2, GET 1 FREE!**

Code: SL3025  
Reg. \$992/ea  
**Net: \$661.33/ea**

**Maxso® Low Speed**

E-Type Connection

**BUY 3, GET 1 FREE!**



1 Year Warranty **Mix and Match**

<b>A. SL2026</b>	1:1 Push Button (CA Burs)	Reg \$168/ea	<b>Net \$126/ea</b>
<b>B. SL2028</b>	1:1 Latch Type (CA Burs)	\$132/ea	<b>\$99/ea</b>
<b>C. SL2034</b>	1:1 Screw-On	\$120/ea	<b>\$90/ea</b>
<b>D. MT2010</b>	Air Motor, 4 Hole, 20,000RPM	\$284/ea	<b>\$213/ea</b>
<b>E. ST2010</b>	1:1 Straight Nose Cone	\$108/ea	<b>\$81/ea</b>
<b>F. SL1002</b>	4:1 Nose Cone, For DPA	\$218/ea	<b>\$163.50/ea</b>

Mix & Match any 3. Promotional goods must be of equal or lesser value than the lowest priced item purchased.

**Slimax-C Plus®**

LED Curing Light System

**BUY 1, GET 1 FREE!**

- Built-In Digital Radiometer
- LED Display Panel
- Cordless, Portable Design
- Cures 850 Times Per Charge
- Up to 2000 mW/cm ± 10%
- Rechargeable Samsung Lithium-Ion Battery
- Solid Aluminum Body
- 1 Year Warranty

Color	Code
Silver	SM1003P-S
Blue	SM1003P-B
Pink	SM1003P-P

Reg. \$628/ea **Net: \$314/ea**

