

AirLight[®] M800 Plus

Generator + LED High Speed

- 23W of Torque
- 410,000 RPM
- Ceramic Ball Bearings
- Titanium Body
- 1 Year Warranty

Backend	Code
4 Hole	HP3008P
Beyes [®]	HP6018P
W&H [®]	HP3022P
Star [®]	HP3054P
NSK [®]	HP3043P
KAVO [®]	HP3036P

Reg. \$548/ea **Net: \$411/ea**

BUY 3, GET 1 FREE

Triple Spray
Direct-LED



Built-in Generator



AirLight[®] M800X

Generator + LED High Speed

- 27W of Torque
- 430,000 RPM
- Instant Stop Capability
- X-Ball Bearings
- Extra Quiet
- Micro-Tex Slip-Proof Coating
- 25 Months Warranty

Backend	Code
Beyes [®]	HP6018X
Star [®]	HP3054X
NSK [®]	HP3043X
KAVO [®]	HP3036X

Reg. \$1,275/ea

Quattro Spray
Direct-LED



Built-in Generator

EZ-Swivel[®]

Couplers

Beyes[®] PD 4 Hole
HP6008
Reg. \$152/ea
Net: \$114/ea



NSK[®] Type 4 Hole, Fiber-optic
HP3046
Reg. \$252/ea **Net: \$189/ea**

NSK[®] Type 6 Hole, Fiber-optic
HP2069
Reg. \$252/ea **Net: \$189/ea**

Kavo[®] Type 4 Hole,
Non-optic
HP2063
Reg. \$252/ea
Net: \$189/ea



Kavo[®] Type 6 Hole,
Fiber-optic
HP2064
Reg. \$252/ea
Net: \$189/ea

Maxso[®] E800

Electric Micromotor System Package

- Powerful: 5.0 Ncm of torque
- Steady Performance:
Speed of rotations is controlled between 2,000 and 40,000min-1
- Smooth Operation:
Torque ripple and heat reduction technology eliminates vibrations
- Convenient & Durable:
Smooth and easy to clean surface with autoclavable brushless motor
- 3 Programmable Memory Settings:
Easily switch between speed presets for more efficient workflow
- New Touch Screen Interface:
Easily and quickly navigate through menu settings on a color LCD display panel
- Warranty: 25 months



MT8001KIT1:
E800P+One X99L Attachment
Reg. \$2,398
Sale: \$2,075



MT8001KIT2:
E800P+Two X99L Attachment
Reg. \$3,118
Sale: \$2,699

Maxso[®] X99L

1:5 Contra Angle

- Anti-Heat System
- Ceramic Ball Bearing

- ISO E-Type Connection
- Quattro Spray
- Fiber-optic
- 1 Year Warranty



Code: SL3025
Reg. \$992/ea
Net: \$661.33/ea

BUY 2, GET 1 FREE

BUY (3) AirLight M800X, GET (1) CareMaster FREE



Handpiece
Maintenance
System

Code: HPL2010

Maxso[®] Low Speed

Internal Water Spray, E-Type Connection

BUY 2, GET 1 FREE

Mix and Match 1 Year Warranty

S20A-IS/E6

- Straight Nose Cone
- 1:1
- Internal Spray
- Fiber-Optic

Code: ST2009
Reg. \$414/ea
Net: \$276/ea



M29L-MH01

- Contra Angle
- Internal Spray
- Fiber-Optic
- Push Button
- CA Burs

Code: SL3068
Reg. \$430/ea
Net: \$286.67/ea



Make it a mix and match with the fine print that promotional goods must be equal or lesser value

Slimax-C Plus[®]

LED Curing Light System

BUY 2, GET 1 FREE

- Powerful up to 2000 mW/cm2 (±10%)
- Built-in Digital Radiometer
- 3 Modes: Full Power, Ramping, Pulsing
- Rechargeable Lithium-ion Battery
- Cordless Design
- 5-Second Curing
- 1 Year Warranty

Color	Code
Silver	SM1003P-S
Blue	SM1003P-B
Pink	SM1003P-PI

Reg. \$660/ea **Net: \$440/ea**



Maxso® M200E

Non-Optic High Speed Handpiece

- 19W of Torque
- 400,000 RPM
- Single Spray
- Push Button
- Low Turbine Noise
- No Swivel
- 6 Months Warranty



Backend Code
4 Hole HP2105

Reg. \$128/ea
Sale: 5 for \$550

Maxso® X200X

Fiber-optic handpiece

- 27W of Torque
- 430,000 RPM
- Instant Stop Capability
- Quattro Spray
- Titanium Body
- 25 Month Warranty



Backend Code
NSK® HP2066X
KAVO® HP2060X

Reg. \$1,048/ea
Net: \$698.67/ea

BUY 2, GET 1 FREE

Maxso® M200 Smart

Non-Optic High Speed Handpiece

- 19W of Torque
- 400,000 RPM
- Push Button Chuck
- Micro-tex Surface
- Triple Water Spray
- 1 Year Warranty



Backend Code
4 Hole HP2009

Reg. \$240/ea
Net: \$180/ea

BUY 3, GET 1 FREE

Maxso® Low Speed

E-Type Connection

BUY 3, GET 1 FREE

Mix and Match 1 Year Warranty



	Reg	Net
A. SL2026 1:1 Push Button (CA Burs)	\$176	\$132
B. SL2028 1:1 Latch Type (CA Burs)	\$147	\$110.25
C. SL2034 1:1 Screw-On	\$137	\$102.75
D. MT2010 Air Motor, 4 Hole, 20,000RPM	\$298	\$223.50
E. ST2010 1:1 Straight Nose Cone	\$114	\$85.50
F. SL1002 4:1 Nose Cone, For DPA	\$278	\$208.50

Mix & Match any 3. Promotional goods must be of equal or lesser value than the lowest priced item purchased.

Maxso® X200

Fiber-optic High Speed Handpiece

- 23W of Torque
- Triple Spray
- Push Button Chuck
- Swivel
- 1 Year Warranty



Backend Code
Kavo® HP2060
NSK® HP2066

Reg. \$480/ea **Net: \$360/ea**

BUY 3, GET 1 FREE

APEXPILOT®

Endo Motor

G2 (AL2060)

- Can be connected to matched apex locator Accurator® A2 (sold separately)
- High Performance Brushless Motor
- 6:1 Super Mini Contra-angle
- Efficient versatile functions



Reg. \$754/ea **Sale: \$597/ea**

G3 (AL2070)

- Built-in Apex Locator
- Integrated Length Determination
- All-angle OLED Screen
- Supports most major brands of NiTi files



Reg. \$1,258/ea **Sale: \$999/ea**

Canaview®

Intraoral Camera

- True Auto-Focus
- Anti-Fog Design
- Responsive Capture Button
- 5 MP High Resolution Images
- Intuitive Camera Lighting System
- No Manual Adjustments Required
- Direct USB Connection to Computer
- Super Light-Weight
- Complimentary Image Management Software



*Compatible with: Abeldent, Adstra, Apixia, Apteryx, CaptureLink, Clio Soft, Cuve Hero, DBSWIN, Dentrix Ascend, Dentrix Image, Dexis, Diagona, Eaglesoft, EVAsoft, ImageFX, Kodak, MiPACS, Mogo, Open Dental, Romexis, Schick, Sidexis, Sopro Imaging, Tiger view, Umbie, VaTech EzDent, VisiQuick, Visix, VixWin

Code: EV2060
Reg. \$904/ea **Net: \$452/ea**



ACCURATOR® A2

Apex Locator, compatible with G2 Endo Motor

More stable, stronger anti-interference ability.
25 National Patents & 8 Invention Patents



(AL2080) Reg. \$838/ea **Sale: \$670/ea**