

## AirLight® M800 X extraordinary

Powerful Quiet Revolutionary

- 27W of Torque
- 430,000 RPM
- Instant Stop Capability
- X-Ball Bearings
- Extra Quiet
- Micro-TEX Slip-Proof Coating
- 25 Months Warranty

Backend	Code
Beyes®	(HP6018X)
Star®	(HP3054X)
NSK®	(HP3043X)
KAVO®	(HP3036X)



Reg: \$1,275/ea Net: \$637<sup>50</sup>/ea

**BUY 1, GET 1 FREE!**

### AirLight® M800 Plus

LED High Speed Handpiece

- 23W of Torque
- 410,000 RPM
- Ceramic Ball Bearings
- Handles up to 60 PSI of air pressure
- Titanium Body
- 1 Year Warranty

Backend	Code
4 Hole	(HP3008P)
Beyes®	(HP6018P)
W&H®	(HP3022P)
Star®	(HP3054P)
NSK®	(HP3043P)
KAVO®	(HP3036P)

**Direct LED Without Fiber Optics!  
Just Brilliantly Bright Illumination!**

Reg. \$548/ea

**Net: \$365<sup>33</sup>/ea**

**BUY 2, GET 1 FREE**



Triple Spray

Direct-LED

Built-in Generator



### Maxso® X200

Fiber-optic High Speed Handpiece

- 23W of Torque
- Triple Spray
- Push Button Chuck
- Swivel
- 1 Year Warranty

Backend	Code
Kavo®	(HP2060)
NSK®	(HP2066)

Reg. \$480/ea

**Net: \$360/ea**

**BUY 3, GET 1 FREE**

Traditional Fiber optics



### Airlight® M800 Surgical 45

LED High Speed Handpiece

- 22W of Torque
- 45° Angled Head
- 410,000 RPM
- Rear Air Exhaust
- Direct LED
- Stream Water
- 1 Year Warranty

Backend	Code
4 Hole	(HP3048)
Beyes®	(HP6038)
NSK®	(HP3044)
KAVO®	(HP3037)
Star®	(HP3055)

Reg. \$598/ea

**Net: \$398<sup>67</sup>/ea**

**BUY 2, GET 1 FREE**

Triple Spray

Direct-LED

Built-in Generator



### EZ-Swivel®

Couplers

Beyes® PD

4 Hole  
Code: (HP6008)

Reg. \$152/ea

**Net: \$101<sup>33</sup>/ea**



NSK® Type

4 Hole  
Code: (HP3046)

6 Hole  
Code: (HP2069)

Reg. \$252/ea

**Net: \$168/ea**



**BUY 2, GET 1 FREE**

Kavo® Type

4 Hole, Non-optic  
Code: (HP2063)

6 Hole, Fiber-optic  
Code: (HP2064)

Reg. \$252/ea

**Net: \$168/ea**



### Slimax-C Plus®

LED Curing Light System

- Built-In Digital Radiometer
- LED Display Panel
- Cordless, Portable Design
- Cures 850 Times Per Charge
- Up to 2000 mW/cm ± 10%
- Rechargeable Samsung Lithium-Ion Battery
- Solid Aluminum Body
- 1 Year Warranty

Color	Code
Silver	(SM1003P-S)
Blue	(SM1003P-B)
Pink	(SM1003P-P)

Reg. \$660/ea

**Net: \$330/ea**



**BUY 1, GET 1 FREE**

### Maxso® M200E

Non-Optic High Speed Handpiece

- 19W of Torque
- 400,000 RPM
- Single Spray
- Push Button
- Low Turbine Noise
- No Swivel
- 6 Months Warranty

Backend	Code
4 Hole	(HP2105)

Reg. \$128/ea



**BUY 5 FOR \$550**

## Maxso® M200 Smart

Non-Optic High Speed Handpiece

- 19W of Torque
- 400,000 RPM
- Push Button Chuck
- Micro-tex Surface
- Triple Water Spray
- 1 Year Warranty

Backend Code  
4 Hole HP2009

Reg. \$240/ea  
**Net: \$180/ea**



**BUY 3, GET 1 FREE**

## Maxso® Air Motor

- Compatible with all standard E-Type Attachments
- Perfectly balanced, lightweight design
- Powerful adjustable speed range 5,000 rpm - 20,000 rpm
- Forward/reverse drive with continuous speed control
- Safe and easy double-lock connection
- Autoclavable
- 1 Year Warranty

Reg. \$298/ea  
**Net: \$198<sup>67</sup>/ea**

Code: (MT2010)



**BUY 2, GET 1 FREE**

## Maxso® Low Speed

E-Type Connection

**BUY 3, GET 1 FREE**

1 Year Warranty



	Reg	Net
A. SL2026 1:1 Push Button (CA Burs) .....	\$176	\$132
B. SL2028 1:1 Latch Type (CA Burs) .....	\$147	\$110.25
C. SL2034 1:1 Screw-On .....	\$137	\$102.75
D. ST2010 1:1 Straight Nose Cone .....	\$114	\$85.50
E. SL1002 4:1 Nose Cone, For DPA .....	\$278	\$208.50

Mix & Match any 3. Promotional goods must be of equal or lesser value than the lowest priced item purchased.

## Canaview®

Intraoral Camera

- True Auto-Focus
- Anti-Fog Design
- Responsive Capture Button
- 5 MP High Resolution Images
- Intuitive Camera Lighting System
- No Manual Adjustments Required
- Direct USB Connection to Computer
- Super Light-Weight
- Complimentary Image Management Software

\*Compatible with: Abeldent, Adstra, Apixia, Apteryx, CaptureLink, Clio Soft, Cuve Hero, DBSWIN, Dentrix Ascend, Dentrix Image, Dexis, Diagona, Eaglesoft, EVAsoft, ImageFX, Kodak, MiPACS, Mogo, Open Dental, Romexis, Schick, Sidexis, Sopro Imaging, Tiger view, Umbie, VaTech EzDent, VisiQuick, Visix, VixWin

Reg. \$904/ea  
**Net: \$452/ea**



Code: (EV2060)

**BUY 1, GET 1 FREE**

## Duray® Intraoral Sensor Package

- Thinnest Sensor 4.4mm
- Durable and Reliable: Waterproof with IP68 protection and can be immersed in disinfectant liquids to prevent cross infections
- High Resolution Imaging: With a theoretical resolution of 25 lp/mm and through advanced CMOS technology, image quality is exceptional
- Compatible with Major Imaging Software: Apeteryx, Dexis, Dentrix Ascend, Eaglesoft, Open Dental, Suni, Schick, Tigerview, Vixwin, XDR, and more!

Model	Code	Reg.	Sale
Size 0	(DS0001)	\$2,220	\$1,999
Size 1	(DS1001)	\$2,940	\$2,450
Size 1.5	(DS1501)	\$3,220	\$2,700
Size 2	(DS2001)	\$3,676	\$3,125
Size 1+2 Package	(DS3001)	\$6,298	\$5,250



24 months warranty



## Duray® V2

DC X-Ray Unit



- Advanced DC Technology
- Excellent Image Quality
- Built-in DAP allows dose monitoring
- Effortless and stable positioning
- 24 Month Warranty

Code: (XR2012)

Reg. \$4,999  
**Sale: \$3,998**



## Duray® Imaging Plate Scanner Package

DC1, Duray Scanner DC1 Package

- Powerful Processing
- Smart Design
- Compact Tabletop Design
- Comfortable
- 24 Months Warranty

Code: (DC1001)

Reg. \$5,550  
**Sale: \$4,999**



**FREE GOODS REDEMPTION:** 1. Go to "http://beyes.ca/en/redeems/new" 2. Fill out the form and submit (Promo code: **USQ424**)

\*\*Free goods must be of equal or lesser value and cannot be combined with any other offers. All offers will require a proof of purchase dated from October 1 - December 31, 2024 and must be reflected on one invoice. Offers must be redeemed no later than January 15, 2025. Please allow 4 to 8 weeks for promotion to be processed. All brands are holders of their respective trademarks.