

## Comfortek® X1

Dental Chair with Delivery Unit

### Includes:

- Stitched upholstery design
- Heavy-duty aluminum backrest
- Double post design
- Flexible delivery system with control panel
- 3 programmable memory positions
- One 6-pin fiber-optic tubing
- Built-in Piezo scaler
- Removable and rotatable stainless steel water outlet
- 4 instrument position on premium swiveling assistant side with premium control panel
- Main switch for all (air, water and electricity)
- 4-way foot switch
- Dentist and assistant stool

Now Upgraded with

- M7 LED Light with Composite Safe Mode
- Extended 8 years warranty

Reg. \$17,599

**Sale: \$14,079**



## AirLight® M800 Smart

LED High Speed Handpiece

- 28W of Torque
- 410,000 RPM
- Ceramic Ball Bearings
- Handles up to 60 PSI of air pressure
- Titanium Body
- 12 Months Warranty

Backend	Code
4 Hole	(HP3008)
Beye®	(HP6018)
W&H®	(HP3022)
Star®	(HP3054)
NSK®	(HP3043)
KAVO®	(HP3036)

Reg. \$812.45/ea

**Net: \$541.<sup>63</sup>/ea**

Triple Spray

Direct-LED

Built-in Generator



## Airlight® M800 Surgical 45

LED High Speed Handpiece

- 25W of Torque
- 45 Degree Angled Head
- 410,000 RPM
- Rear Air Exhaust
- Direct LED
- Stream Water
- 12 Months Warranty

Backend	Code
4 Hole	(HP3048)
Beyes®	(HP6038)
NSK®	(HP3044)
KAVO®	(HP3037)
Star®	(HP3055)

Reg. \$883.95/ea

**Net: \$589.<sup>30</sup>/ea**

Triple Spray

Direct-LED

Built-in Generator



## AirLight® M800X

LED High Speed Handpiece

- 35W of Torque
- 430,000 RPM
- Extra Quiet
- 25 Months Warranty

**Torque 35w**

Backend	Code
Beyes®	(HP6018X)
Star®	(HP3054X)
NSK®	(HP3043X)
KAVO®	(HP3036X)

Reg. \$1,884.95/ea **Net: \$942.<sup>48</sup>/ea**

Limitless PSI Technology

Experience the freedom from operating pressure control

X Ball Bearings  
Lasts 4x longer than other ball bearings

Instant Stop Technology

Eliminates backflow and prevents cross-contamination



## Slimax-C Plus®

LED Curing Light System

- Built-In Digital Radiometer
- LED Display Panel
- Cordless, Portable Design
- Cures 850 Times Per Charge
- Up to 2000 mW/cm ± 10%
- Rechargeable Samsung Lithium-Ion Battery
- Solid Aluminum Body
- 12 Months Warranty

Ultra-fast  
5-Second  
Cure Time!

Color	Code
Silver	(SM1003P-S)
Blue	(SM1003P-B)
Pink	(SM1003P-P)

Reg. \$812.45/ea **Net: \$406.<sup>23</sup>/ea**



## Maxso® X200X

Fiber-Optic High Speed Handpiece

- 35W of Torque
- 430,000 RPM
- Instant Stop Capability
- Quattro Spray
- Titanium Body
- 25 Month Warranty

**Torque 35w**

Backend	Code
NSK®	HP2066X
KAVO®	HP2060X

Reg. \$1,559.95/ea **Net: \$779.<sup>99</sup>/ea**



## Canaview®

Intraoral Camera

- 5MP High resolution images
- True auto-focus and anti-fog lens
- Plug-and-play with direct USB connection
- Intuitive camera lighting system
- Super lightweight
- Compatible with all major imaging software including Dentrax, Dexis, Eaglesoft, Curve Hero, Open Dental ...and more!
- 12 Months Warranty

Reg. \$1,298.65/ea **Net: \$865.<sup>77</sup>/ea** Code: (EV2060)



**FREE GOODS REDEMPTION: 1.** Go to “<http://beyesdental.com/en/redeems/new>” **2.** Fill out the form and submit (Promo code: **CAQ126**)

\*Free goods must be of equal or lesser value and cannot be combined with any other offers. All offers will require a proof of purchase dated from January 1 - March 31, 2026 and must be reflected on one invoice. Offers must be redeemed no later than April 15, 2026. Please allow 4 to 8 weeks for promotion to be processed. All brands are holders of their respective trademarks. All promotions are subject to change or cancellation without notice. Beyes and its authorized distributors reserve the right to modify offers, correct errors and limit quantities at any time; in the event of a discrepancy, the price and terms confirmed at the point of sale will prevail.

## Maxso® E800

Electric Handpiece System

- **Powerful:** 5.0 Ncm of torque
- **Steady Performance:**  
Speed of rotations is controlled between 2,000 and 40,000min-1
- **Smooth Operation:**  
Torque ripple and heat reduction technology eliminates vibrations
- **Convenient & Durable:**  
Smooth and easy to clean surface with autoclavable brushless motor
- **3 Programmable Memory Settings:**  
Easily switch between speed presets for more efficient workflow
- **New Touch Screen Interface:**  
Easily and quickly navigate through menu settings on a color LCD display panel
- **Warranty:** 25 Months

## REVOLUTIONARY TECHNOLOGY MADE EASY!



**Code:** MT8001KIT1  
**(1) E800P + (1) X99L**

Reg. \$3,368.75  
**Sale: \$2,832.<sup>90</sup>**



**Code:** MT8001KIT2  
**(1) E800P + (2) X99L**

Reg. \$4,370.65  
**Sale: \$3,684.<sup>29</sup>**

## Maxso® X99L

1:5 Contra Angle

- Anti-Heat System
- Quattro Spray
- Ceramic Ball Bearing
- Fiber-optic
- ISO E-Type Connection
- 12 Months Warranty

**Code:** SL3025

Reg. \$1,416.95/ea  
**Net: \$944.<sup>63</sup>/ea**



**BUY 2, GET 1 FREE**

## EZ-Swivel®

Couplers

Mix & Match

**BUY 2, GET 1 FREE**

**Beyes® PD**  
**4 Hole**  
**Code:** (HP6008)  
Reg. \$220.95/ea  
**Net: \$147.<sup>30</sup>/ea**



**NSK® Type**  
**4 Hole**  
**Code:** (HP3046)  
**6 Hole**  
**Code:** (HP2069)  
Reg. \$370.45/ea  
**Net: \$246.<sup>97</sup>/ea**



**Kavo® Type**  
**4 Hole, Non-optic**  
**Code:** (HP2063)  
**6 Hole, Fiber-optic**  
**Code:** (HP2064)  
Reg. \$370.45/ea  
**Net: \$246.<sup>97</sup>/ea**



Mix & Match any 2. Promotional goods must be of equal or lesser value than the lowest priced item purchased.

## Duray® Intraoral Sensor Package



- Thinnest Sensor 4.4mm
- Durable and Reliable: Waterproof with IP68 protection and can be immersed in disinfectant liquids to prevent cross infections
- High Resolution Imaging: With a theoretical resolution of 25 lp/mm and through advanced CMOS technology, image quality is exceptional
- Compatible with Major Imaging Software: Apeteryx, Dexis, Dentrax Ascend, Eaglesoft, Open Dental, Suni, Schick, Tigerview, Vixwin, XDR, and more!
- 24 Months Warranty

Model	Code	Reg.	Sale
Size 0	(DS0001)	\$3,118.65	<b>\$2,865.18</b>
Size 1	(DS1001)	\$4,159.95	<b>\$3,617.26</b>
Size 1.5	(DS1501)	\$4,549.95	<b>\$3,822.00</b>
Size 2	(DS2001)	\$5,198.65	<b>\$4,436.29</b>
Size 1+2 Package	(DS3001)	\$8,969.95	<b>\$7,644.00</b>



## Maxso® Low Speed

E-Type Connection

**BUY 3, GET 1 FREE**

Mix & Match



12 Months Warranty

**Reg**

**Net**

<b>A. SL2026</b>	1:1 Push Button (CA Burs).....	\$246.95	<b>\$185.21</b>
<b>B. SL2028</b>	1:1 Latch Type (CA Burs).....	\$214.45	<b>\$160.84</b>
<b>C. SL2034</b>	1:1 Screw-On .....	\$194.95	<b>\$146.21</b>
<b>D. ST2010</b>	1:1 Straight Nose Cone.....	\$188.45	<b>\$141.34</b>
<b>E. SL1002</b>	4:1 Nose Cone, For DPA .....	\$402.95	<b>\$302.21</b>

Mix & Match any 3. Promotional goods must be of equal or lesser value than the lowest priced item purchased.

## CareMaster®

Automatic Handpiece Maintenance System

Maintain Your Most Significant Investments in the Office!  
Clean 3 Handpieces Per Cycle!

- Maximize the lifespan of handpieces
- Economical, efficient and time-saving
- Filtered air used to rotate the handpieces
- Detachable door for easy cleaning
- Provides care for all brands of handpieces
- 12 Months Warranty



\* Requires Airline Connection  
\* Coupler Sold Separately

**Code:** HPL2010  
2 Turbines & 1 E-Type Attachment


Reg. \$1,754.95/ea  
**Sale: \$1,559.<sup>80</sup>/ea**

**Code:** HPL2011  
1 Turbines & 2 E-Type Attachments

Reg. \$1,917.45/ea  
**Sale: \$1,754.<sup>30</sup>/ea**

## Duray® 2D, Ceph & 3D

Extraoral X-ray Systems



### Advanced Imaging for Every Practice

- Quick Cephalometrics with Touchscreen Access
- 2D, Ceph, and Fully Upgradeable 3D Options
- Automatic Metal Artifact & Noise Reduction (AMAR & ANR)
- Fast & Easy Patient Positioning | Plug & Play Installation
- Smart Backlight | Image Mirroring | Advanced Sensor
- 5 Year Warranty Across All Models

Model	Code
Duray 2D Pan	(XR3020)
Duray 2D Ceph	(XR3021)
Duray 3D	(XR3030)
Duray 3D Ceph	(XR3031)



**Contact us for  
complete specs  
and pricing**

www.beyesdental.com  
info@beyes.ca  
1-855-603-1888

**FREE GOODS REDEMPTION: 1.** Go to "http://beyesdental.com/en/redeems/new" **2.** Fill out the form and submit (Promo code: **CAQ126**)

\*Free goods must be of equal or lesser value and cannot be combined with any other offers. All offers will require a proof of purchase dated from January 1 - March 31, 2026 and must be reflected on one invoice. Offers must be redeemed no later than April 15, 2026. Please allow 4 to 8 weeks for promotion to be processed. All brands are holders of their respective trademarks. All promotions are subject to change or cancellation without notice. Beyes and its authorized distributors reserve the right to modify offers, correct errors and limit quantities at any time; in the event of a discrepancy, the price and terms confirmed at the point of sale will prevail.