

Canaview® Intraoral Camera

- True Auto-Focus
- Anti-Fog Design
- Responsive Capture Button
- 5 MP High Resolution Images
- Direct USB Connection to Computer
- Intuitive Camera Lighting System
- Super Lightweight
- Complimentary Image Management Software
- 12 Months Warranty

***Compatible with:** Abeldent, Adstra, Apixia, Apteryx, CaptureLink, Clio Soft, Cuve Hero, DBSWIN, Dentrax Ascend, Dentrax Image, Dexis, Diagora, Eaglesoft, EVAsoft, ImageFX, Kodak, MiPACS, Mogo, Open Dental, Romexis, Schick, Sidexis, Sopro Imaging, Tiger view, Umbie, VaTech EzDent, VisiQuick, Visix, VixWin...and more!



Code: (EV2060)

BUY 1, GET 1 FREE

Reg. \$1,287/ea **Net: \$643.50/ea**

AirLight® M800X

LED High Speed Handpiece

- 30W of Torque
- 430,000 RPM
- Extra Quiet
- 25 Months Warranty

| Backend | Code |
|---------|-----------|
| Beyes® | (HP6018X) |
| Star® | (HP3054X) |
| NSK® | (HP3043X) |
| KAVO® | (HP3036X) |

Reg. \$1,818.70/ea

Net: \$1,212.47/ea

**Torque
30W**

Limitless PSI Technology

Experience the freedom from operating pressure control

X Ball Bearings

Lasts 4x longer than other ball bearings

Instant Stop Technology

Eliminates backflow and prevents cross-contamination



BUY 2, GET 1 FREE

Airlight® M800 Surgical 45

LED High Speed Handpiece

- 22W of Torque
- 45 Degree Angled Head
- 410,000 RPM
- Rear Air Exhaust
- Direct LED
- Stream Water
- 12 Months Warranty

| Backend | Code |
|---------|----------|
| 4 Hole | (HP3048) |
| Beyes® | (HP6038) |
| NSK® | (HP3044) |
| KAVO® | (HP3037) |
| Star® | (HP3055) |

Reg. \$845/ea

Net: \$563.33/ea

**Triple Spray
Direct-LED**



**Built-in
Generator**

AirLight® M800 Smart

LED High Speed Handpiece

- 23W of Torque
- 410,000 RPM
- Ceramic Ball Bearings
- Handles up to 60 PSI of air pressure
- Titanium Body
- 12 Months Warranty

| Backend | Code |
|---------|--------|
| 4 Hole | HP3008 |
| Beye® | HP6018 |
| W&H® | HP3022 |
| Star® | HP3054 |
| NSK® | HP3043 |
| KAVO® | HP3036 |

Reg. \$778.70/ea

Net: \$584.03/ea

Triple Spray

Direct-LED

**Built-in
Generator**



BUY 3, GET 1 FREE

EZ-Swivel®

Couplers



Beyes® PD

4 Hole
Code: HP6008
Reg. \$214.50/ea

Net: \$160.88/ea



Kavo® Type

4 Hole, Non-optic
Code: HP2063
Reg. \$357.50/ea

Net: \$268.13/ea



NSK® Type

4 Hole, Fiber-optic
Code: HP3046
Reg. \$357.50/ea

Net: \$268.13/ea



Kavo® Type

6 Hole, Fiber-optic
Code: HP2064
Reg. \$357.50/ea

Net: \$268.13/ea

BUY 3, GET 1 FREE

Mix & Match

Mix & Match any 3. Promotional goods must be of equal or lesser value than the lowest priced item purchased.

Maxso® X200X

Fiber-Optic High Speed Handpiece

- 27W of Torque
- 430,000 RPM
- Instant Stop Capability
- Quattro Spray
- Titanium Body
- 25 Months Warranty

| Backend | Code |
|---------|---------|
| NSK® | HP2066X |
| KAVO® | HP2060X |

Reg. \$1,495/ea

Net: \$747.50/ea



BUY 1, GET 1 FREE

Slimax-C Plus®

LED Curing Light System

- Built-In Digital Radiometer
- LED Display Panel
- Cordless, Portable Design
- Cures 850 Times Per Charge
- Up to 2000 mW/cm ± 10%
- Rechargeable Samsung Lithium-Ion Battery
- Solid Aluminum Body
- 12 Months Warranty

| Color | Code |
|--------|--------------|
| Silver | (SM1003P-S) |
| Blue | (SM1003P-B) |
| Pink | (SM1003P-PI) |

Reg. \$908.70/ea

Net: \$605.80/ea



BUY 2, GET 1 FREE

FREE GOODS REDEMPTION: 1.Go to "http://beyesdental.com/en/redeems/new" 2.Fill out the form and submit (Promo code: CAQ325)

*Free goods must be of equal or lesser value and cannot be combined with any other offers. All offers will require a proof of purchase dated from July 1 - September 30, 2025 and must be reflected on one invoice. Offers must be redeemed no later than October 15, 2025. Please allow 4 to 8 weeks for promotion to be processed. All brands are holders of their respective trademarks.

Maxso® M200E

Non-Optic High Speed Handpiece

- 19W of Torque
- 400,000 RPM
- Single Spray
- Push Button
- Low Turbine Noise
- No Swivel
- 6 Months Warranty

Backend 4 Hole
Code HP2105

Reg. \$205.40/ea
Net: \$171⁻¹⁷/ea

BUY 5, GET 1 FREE



Maxso® M200 Smart

Non-Optic High Speed Handpiece

- 19W of Torque
- 400,000 RPM
- Push Button Chuck
- Micro-tex Surface
- Triple Water Spray
- 12 Months Warranty

Backend 4 Hole
Code HP2009

Reg. \$336.70/ea
Net: \$252⁻⁵³/ea

BUY 3, GET 1 FREE



Maxso® Low Speed

E-Type Connection

BUY 3, GET 1 FREE

Mix & Match



12 Months Warranty

| | | Reg | Net |
|------------------|---------------------------------|----------|-----------------|
| A. SL2026 | 1:1 Push Button (CA Burs) | \$247 | \$185.25 |
| B. SL2028 | 1:1 Latch Type (CA Burs) | \$208 | \$156 |
| C. SL2034 | 1:1 Screw-On | \$195 | \$146.25 |
| D. ST2010 | 1:1 Straight Nose Cone | \$180.70 | \$135.53 |
| E. SL1002 | 4:1 Nose Cone, For DPA | \$388.70 | \$291.53 |
| F. MT2010 | Air Motor, 4 Holes | \$422.50 | \$316.88 |

Mix & Match any 3. Promotional goods must be of equal or lesser value than the lowest priced item purchased.

Duray® Intraoral Sensor Package

- Thinnest Sensor 4.4mm
- Durable and Reliable: Waterproof with IP68 protection and can be immersed in disinfectant liquids to prevent cross infections
- High Resolution Imaging: With a theoretical resolution of 25 lp/mm and through advanced CMOS technology, image quality is exceptional
- Compatible with Major Imaging Software: Apeteryx, Dexis, Dentrix Ascend, Eaglesoft, Open Dental, Suni, Schick, Tigerview, Vixwin, XDR, and more!
- 24 Months Warranty



| Model | Code | Reg. | Sale |
|------------------|----------|------------|-------------------|
| Size 0 | (DS0001) | \$3,003 | \$2,728.70 |
| Size 1 | (DS1001) | \$4,011.80 | \$3,445 |
| Size 1.5 | (DS1501) | \$4,394 | \$3,640 |
| Size 2 | (DS2001) | \$5,005 | \$4,225 |
| Size 1+2 Package | (DS3001) | \$8,580 | \$7,280 |



Duray® 2D & Ceph

Extraoral X-ray Series

- Smart Backlight
- Integrated Touchscreen
- Advanced Sensor
- New Dynamic and Automatic Image Focusing
- One Shot Ceph
- Advanced Sensor
- Smart Backlight

Duray 2D Pan
Code: XR3020

Duray 2D Ceph
Code: XR3021




Duray® 3D & Ceph

Extraoral X-ray Series

- Touchscreen Next To Patient
- Quick Cephalometrics
- Instantly Upgradeable
- Automatic Metal Artifact Reduction (AMAR)
- Automatic Noise Reduction (ANR)
- Advanced Technology
- Mirroring
- Smart Backlight
- Safe & Easy
- Plug & Play
- Quick Positioning

Duray 3D Pan
Code: XR3030

Duray 3D Ceph
Code: XR3031



Accurator® A2

Apex Locator, Compatible with G2 Endo Motor

- 4.5" wide, color LCD screen provides clear images with high contrast
- Easy-to-use interface with adjustable inclination
- Accurate measurement in wet and dry canals
- Automatic calibration
- Autoclavable file clip, file hook, file probe
- 12 Months Warranty

Code: AL2080

Reg. \$1,198.34 **Sale: \$955⁻⁵⁰**



Duray® V2

DC X-Ray Unit

- Advanced DC Technology
- Excellent Image Quality
- Built-in DAP allows dose monitoring
- Effortless and stable positioning
- 24 Months Warranty

Code: (XR2012)

Reg. \$5,717.40 **Sale: \$5,193⁻⁵⁰**



FREE GOODS REDEMPTION: 1.Go to "http://beyes.ca/en/redeems/new" 2.Fill out the form and submit (Promo code: **CAQ325**)

*Free goods must be of equal or lesser value and cannot be combined with any other offers. All offers will require a proof of purchase dated from July 1 - September 30, 2025 and must be reflected on one invoice. Offers must be redeemed no later than October 15, 2025. Please allow 4 to 8 weeks for promotion to be processed. All brands are holders of their respective trademarks.