

# AirLight<sup>®</sup> M800 X extraordinary

Powerful Quiet Revolutionary

- 27W of Torque
- 430,000 RPM
- Instant Stop Capability
- X-Ball Bearings
- Extra Quiet
- Micro-Tex Slip-Proof Coating
- 25 Months Warranty

Backend	Code
Beyes <sup>®</sup>	(HP6018X)
Star <sup>®</sup>	(HP3054X)
NSK <sup>®</sup>	(HP3043X)
KAVO <sup>®</sup>	(HP3036X)



Reg: \$1,498/ea **Net: \$749/ea**

**BUY 1, GET 1 FREE!**

## AirLight<sup>®</sup> M800 Smart

LED High Speed Handpiece

- 23W of Torque
- 410,000 RPM
- Ceramic Ball Bearings
- Handles up to 60 PSI of air pressure
- Titanium Body
- 1 Year Warranty

Backend	Code
4 Hole Beyes <sup>®</sup>	(HP3008)
W&H <sup>®</sup>	(HP3022)
Star <sup>®</sup>	(HP3054)
NSK <sup>®</sup>	(HP3043)
KAVO <sup>®</sup>	(HP3036)

**Direct LED Without Fiber Optics!  
Just Brilliantly Bright Illumination!**

Reg. \$678<sup>60</sup>/ea  
**Net: \$452<sup>40</sup>/ea**

**BUY 2, GET 1 FREE**



## Maxso<sup>®</sup> X200

Fiber-optic High Speed Handpiece

- 23W of Torque
- Triple Spray
- Push Button Chuck
- Swivel
- 1 Year Warranty

Backend	Code
Kavo <sup>®</sup>	(HP2060)
NSK <sup>®</sup>	(HP2066)

Reg. \$626/ea **Net: \$469<sup>50</sup>/ea**

**BUY 3, GET 1 FREE**

Traditional Fiber optics



## AirLight<sup>®</sup> M800 Surgical 45

LED High Speed Handpiece

- 22W of Torque
- 45 Degree Angled Head
- 410,000 RPM
- Rear Air Exhaust
- Direct LED
- Stream Water
- 1 Year Warranty

Backend	Code
4 Hole Beyes <sup>®</sup>	(HP3048)
NSK <sup>®</sup>	(HP3044)
KAVO <sup>®</sup>	(HP3037)
Star <sup>®</sup>	(HP3055)

Reg. \$796/ea  
**Net: \$530<sup>67</sup>/ea**



**BUY 2, GET 1 FREE**

## EZ-Swivel<sup>®</sup>

Couplers Mix & Match

\*of the same or lesser value of cheapest item purchased.

Beyes <sup>®</sup> PD	Beyes <sup>®</sup> QD	NSK <sup>®</sup> Type	Kavo <sup>®</sup> Type
B2, Non-optic (HP6002)	B2, Non-optic (HP3019)	4 Hole, Non-optic (HP3046)	4 Hole, Non-optic (HP2063)
4 Hole, Non-optic (HP6008)	4 Hole, Non-optic (HP3018)	6 Hole, Fiber-optic (HP2069)	6 Hole, Fiber-optic (HP2064)
Reg. \$186/ea <b>Net: \$139.50/ea</b>	Reg. \$186/ea <b>Net: \$139.50/ea</b>	Reg. \$328/ea <b>Net: \$246/ea</b>	Reg. \$328/ea <b>Net: \$246/ea</b>

**BUY 3, GET 1 FREE\***

## Slimax-C Plus<sup>®</sup>

LED Curing Light System

- Built-In Digital Radiometer
- LED Display Panel
- Cordless, Portable Design
- Cures 850 Times Per Charge
- Up to 2000 mW/cm ± 10%
- Rechargeable Samsung Lithium-Ion Battery
- Solid Aluminum Body
- 1 Year Warranty

Color	Code	Reg.	Net:
Silver	(SM1003P-S)	\$838/ea	<b>\$419/ea</b>
Blue	(SM1003P-B)		
Pink	(SM1003P-P)		



## Duray<sup>®</sup> Intraoral Sensor Package



- Thinnest Sensor 4.4mm
- Durable and Reliable: Waterproof with IP68 protection and can be immersed in disinfectant liquids to prevent cross infections
- High Resolution Imaging: With a theoretical resolution of 25 lp/mm and through advanced CMOS technology, image quality is exceptional
- Compatible with Major Imaging Software: Apeteryx, Dexas, Dentrix Ascend, Eaglesoft, Open Dental, Suni, Schick, Tigerview, Vixwin, XDR, and more!

Model	Code	Reg.	Sale
Size 0	(DS0001)	\$3,036	<b>\$2,599</b>
Size 1	(DS1001)	\$3,900	<b>\$3,185</b>
Size 1.5	(DS1501)	\$4,200	<b>\$3,510</b>
Size 2	(DS2001)	\$4,790	<b>\$4,065</b>
Sensor Package (Size 1+Size2)	(DS3001)	\$8,190	<b>\$6,825</b>

24 months warranty



**FREE GOODS REDEMPTION:** 1. Go to "http://beyes.ca/en/redeems/new" 2. Fill out the form and submit (Promo code: CAQ424)

\*\*Free goods must be of equal or lesser value and cannot be combined with any other offers. All offers will require a proof of purchase dated from October 1 to December 31, 2024 and must be reflected on one invoice. Offers must be redeemed no later than January 15, 2025. Please allow 4 to 8 weeks for promotion to be processed. All brands are holders of their respective trademarks.

**CI5-E1**


Chemical Integrator for Steam Sterilization

- Type 5 Test, Single-use moving front chemical indicators for monitoring Steam sterilization process.
- Consists of a special wick and a steam and temperature sensitive chemical pellet.
- The migration is visible through a zone marked Accept or Reject.

**4+1 FREE**

Reg. \$29/pack

**Net: \$23<sup>.20</sup>/pack**



**Code: CI1200**  
100/Pack



**CI5-P1**

Unique Point Chemical Integrator for Steam Sterilization

- Type 5 Test, Single-use unique point chemical indicators that consist of paper strips printed with indicator ink.
- Intended to be placed in each pack, pouch, container, tray or other containment device to function as an independent monitor of critical parameters for sterilization cycles.
- Colour change indicates effective sterilization processes.

**4+1 FREE**

Reg. \$29/pack

**Net: \$23<sup>.20</sup>/pack**



**Code: CI1100**  
200/Pack



**CI2-P1**

Bowie-Dick Test Card

- Bowie Dick Test, Consists of a Type 2, metals free chemical indicator. Used in Class B Sterilizers.
- Use daily to monitor the effectiveness of air removal.
- Cycle information can be recorded on the back of the indicator.

**4+1 FREE**

Reg. \$239/pack

**Net: \$191<sup>.20</sup>/pack**



**Code: CI1300**  
120/Pack



**Bowie-Dick Test Card Holder**

- For Bowie Dick Test Card
- Metal Card Holder

**4+1 FREE**

Reg. \$110/pack

**Net: \$88/pack**

**Code: CI1400**



**Sensibles**<sup>™</sup>

Universal Sensor Holder Starter Kit

Reg. \$270.40 **Sale: \$205<sup>.20</sup>**

- Everything you need to capture the perfect image is inside!
- Adjusts easily; holds securely
- Perfect for Anterior, Posterior and Bite-Wing images
- Autoclavable and reusable (check manual for details and instructions)
- Durable and 100% latex-free

**COMPATIBLE WITH ALL BRAND SENSORS!!**




**Code: DS1103**

**Include:**

- Positioning Arm Posterior x1
- Positioning Arm Bite Wing x1, Positioning Arm Anterior x1
- Bite Block Large Size x2, Bite Block Medium Size x2
- Posterior Aiming Ring x1, Anterior/Bite Wing Aiming Ring x1

Using Sensibles is as easy as 1-2-3



**1** Attach Positioning Bar



**2** Add Aiming Ring



**3** Snap In Sensor



For Bite-Wings, adjust vertical position