

Comfortek® X1

Dental Chair with Delivery Unit

Includes:

- Stitched upholstery design
- Heavy-duty aluminum backrest
- Double post design
- Flexible delivery system with control panel
- 3 programmable memory positions
- Built-in Piezo scaler
- Removable and rotatable stainless steel water outlet
- 4 instrument position on premium swiveling assistant side with premium control panel
- Main switch for all (air, water and electricity)
- 4-way foot switch
- Dentist and assistant stool

Now Upgraded with

- M7 LED Light with Composite Safe Mode
- Extended 8 years warranty

Reg. \$17,599

Sale: \$14,079



AirLight® M800 Smart

LED High Speed Handpiece

- 23W of Torque
- 410,000 RPM
- Ceramic Ball Bearings
- Handles up to 60 PSI of air pressure
- Titanium Body
- 12 Months Warranty

Backend	Code
4 Hole	(HP3008)
Beye®	(HP6018)
W&H®	(HP3022)
Star®	(HP3054)
NSK®	(HP3043)
KAVO®	(HP3036)

Reg. \$778.70/ea

Net: \$519⁻¹⁵/ea

Triple Spray

Direct-LED

Built-in Generator



BUY 2, GET 1 FREE

Airlight® M800 Surgical 45

LED High Speed Handpiece

- 22W of Torque
- 45 Degree Angled Head
- 410,000 RPM
- Rear Air Exhaust
- Direct LED
- Stream Water
- 12 Months Warranty

Backend	Code
4 Hole	(HP3048)
Beyes®	(HP6038)
NSK®	(HP3044)
KAVO®	(HP3037)
Star®	(HP3055)

Reg. \$845/ea

Net: \$563⁻³³/ea

Triple Spray

Direct-LED

Built-in Generator



BUY 2, GET 1 FREE

AirLight® M800X

LED High Speed Handpiece

- 30W of Torque
- 430,000 RPM
- Extra Quiet
- 25 Months Warranty

Backend	Code
Beyes®	(HP6018X)
Star®	(HP3054X)
NSK®	(HP3043X)
KAVO®	(HP3036X)

Reg. \$1818.70/ea

Net: \$909⁻³⁵/ea

Limitless PSI Technology

Experience the freedom from operating pressure control

X Ball Bearings

Lasts 4x longer than other ball bearings

Instant Stop Technology

Eliminates backflow and prevents cross-contamination



BUY 1, GET 1 FREE

Slimax-C Plus®

LED Curing Light System

- Built-In Digital Radiometer
- LED Display Panel
- Cordless, Portable Design
- Cures 850 Times Per Charge
- Up to 2000 mW/cm ± 10%
- Rechargeable Samsung Lithium-Ion Battery
- Solid Aluminum Body
- 12 Months Warranty

Reg. \$908.70/ea

Net: \$454⁻³⁵/ea

Ultra-fast 5-Second Cure Time!

Color
Silver
Blue
Pink

Code
(SM1003P-S)
(SM1003P-B)
(SM1003P-P)



BUY 1, GET 1 FREE

Maxso® M200E

Non-Optic High Speed Handpiece

- 19W of Torque
- 400,000 RPM
- Single Spray
- Push Button
- Low Turbine Noise
- No Swivel
- 6 Months Warranty

Backend	Code
4 Hole	(HP2105)

Reg. \$205.40/ea

Net: \$171⁻¹⁷/ea

RATED #1



BUY 5, GET 1 FREE

Canaview®

Intraoral Camera

- 5MP High resolution images
- True auto-focus and anti-fog lens
- Plug-and-play with direct USB connection
- Intuitive camera lighting system
- Super lightweight
- Compatible with all major imaging software including Dentrax, Dexis, Eaglesoft, Curve Hero, Open Dental ...and more!
- 12 Months Warranty

Reg. \$1,287/ea

Net: \$643⁻⁵⁰/ea

Code: (EV2060)



FREE GOODS REDEMPTION: 1. Go to "http://beyesdental.com/en/redeems/new" **2.** Fill out the form and submit (Promo code: **CAQ425**)

*Free goods must be of equal or lesser value and cannot be combined with any other offers. All offers will require a proof of purchase dated from October 1 - December 31, 2025 and must be reflected on one invoice. Offers must be redeemed no later than January 15, 2026. Please allow 4 to 8 weeks for promotion to be processed. All brands are holders of their respective trademarks.

Maxso® M200 Smart

Non-Optic High Speed Handpiece

- 19W of Torque
- 400,000 RPM
- Push Button Chuck
- Micro-tex Surface
- Triple Water Spray
- 12 Months Warranty

Backend Code
4 Hole (HP2009)

Reg. \$336.70/ea
Net: \$252^{.53}/ea



BUY 3, GET 1 FREE

Maxso® X200 Smart

Fiber-Optic High Speed Handpiece

- Push Button Chuck
- Micro-tex Surface
- Triple Spray
- Connect directly to Kavo® Multiflex Coupler
- 12 Months Warranty

Backend Code
NSK® (HP2066)
KAVO® (HP2060)

Reg. \$682.50/ea
Net: \$511^{.88}/ea



BUY 3, GET 1 FREE

EZ-Swivel®

Couplers

Beyes® PD
4 Hole
Code: (HP6008)
Reg. \$214.50/ea
Net: \$143/ea



NSK® Type
4 Hole
Code: (HP3046)
6 Hole
Code: (HP2069)
Reg. \$357.50/ea
Net: \$238^{.33}/ea



Kavo® Type
4 Hole, Non-optic
Code: (HP2063)
6 Hole, Fiber-optic
Code: (HP2064)
Reg. \$357.50/ea
Net: \$238^{.33}/ea



Mix & Match any 3. Promotional goods must be of equal or lesser value than the lowest priced item purchased.

Mix & Match

BUY 2, GET 1 FREE

Maxso® X99L

1:5 Contra Angle

- Anti-Heat System
- Ceramic Ball Bearing
- ISO E-Type Connection
- Quattro Spray
- Fiber-optic
- 12 Months Warranty

Code: (SL3025)
Reg. \$1,417/ea
Net: \$944^{.67}/ea



BUY 2, GET 1 FREE

Maxso® Low Speed

E-Type Connection

Mix & Match

A



B



C



D



E



12 Months Warranty

		Reg	Net
A. SL2026	1:1 Push Button (CA Burs).....	\$247	\$185.25
B. SL2028	1:1 Latch Type (CA Burs).....	\$208	\$156
C. SL2034	1:1 Screw-On	\$195	\$146.25
D. ST2010	1:1 Straight Nose Cone.....	\$180.70	\$135.53
E. SL1002	4:1 Nose Cone, For DPA.....	\$388.70	\$291.53

Mix & Match any 3. Promotional goods must be of equal or lesser value than the lowest priced item purchased.

Maxso® E800

Electric Micromotor System Package



MT8001KIT1:
E800P + One X99L Attachment

Reg. \$3,368.82
Sale: \$2,358^{.17}



MT8001KIT2:
E800P + Two X99L Attachment

Reg. \$4,758
Sale: \$3,357^{.17}

Maxso® Air Motor

BUY 2, GET 1 FREE

- Compatible with all standard E-Type Attachments
- Perfectly balanced, lightweight design
- Powerful adjustable speed range 5,000 rpm - 20,000 rpm
- Fast & Easy Patient Positioning | Plug & Play Installation
- Smart Backlight | Image Mirroring | Advanced Sensor
- 5 Year Warranty Across All Models

Reg. \$422.50/ea
Net: \$281^{.67}/ea

Code: (MT2010)



Duray® Intraoral Sensor Package

- Thinnest Sensor 4.4mm
- Durable and Reliable: Waterproof with IP68 protection and can be immersed in disinfectant liquids to prevent cross infections
- High Resolution Imaging: With a theoretical resolution of 25 lp/mm and through advanced CMOS technology, image quality is exceptional
- Compatible with Major Imaging Software: Apertyx, Dexas, Dentrix Ascend, Eaglesoft, Open Dental, Suni, Schick, Tigerview, Vixwin, XDR, and more!
- 24 Months Warranty



Model	Code	Reg.	Sale
Size 0	(DS0001)	\$3,003	\$2,728.70
Size 1	(DS1001)	\$4,011.80	\$3,445
Size 1.5	(DS1501)	\$4,394	\$3,640
Size 2	(DS2001)	\$5,005	\$4,225
Size 1+2 Package	(DS3001)	\$8,580	\$7,280



Duray® 2D, Ceph & 3D

Extraoral X-ray Systems


Advanced Imaging for Every Practice

- Quick Cephalometrics with Touchscreen Access
- 2D, Ceph, and Fully Upgradeable 3D Options
- Automatic Metal Artifact & Noise Reduction (AMAR & ANR)
- Fast & Easy Patient Positioning | Plug & Play Installation
- Smart Backlight | Image Mirroring | Advanced Sensor
- 5 Year Warranty Across All Models

Model	Code
Duray 2D Pan	(XR3020)
Duray 2D Ceph	(XR3021)
Duray 3D	(XR3030)
Duray 3D Ceph	(XR3031)

Contact us for complete specs and pricing

www.beyesdental.com
info@beyes.ca
1-855-603-1888



FREE GOODS REDEMPTION: 1.Go to "http://beyesdental.com/en/redeems/new" 2.Fill out the form and submit (Promo code: **CAQ425**)

*Free goods must be of equal or lesser value and cannot be combined with any other offers. All offers will require a proof of purchase dated from October 1 - December 31, 2025 and must be reflected on one invoice. Offers must be redeemed no later than January 15, 2026. Please allow 4 to 8 weeks for promotion to be processed. All brands are holders of their respective trademarks.

CI5-E1

Chemical Integrator for Steam Sterilization

- Type 5 Test, Single-use moving front chemical indicators for monitoring Steam sterilization process.
- Consists of a special wick and a steam and temperature sensitive chemical pellet.
- The migration is visible through a zone marked Accept or Reject.

100/Pack

Code: (CI1200)

Reg. \$31

Sale: \$25



CI5-P1

Unique Point Chemical Integrator for Steam Sterilization


- Type 5 Test, Single-use unique point chemical indicators that consist of paper strips printed with indicator ink.
- Intended to be placed in each pack, pouch, container, tray or other containment device to function as an independent monitor of critical parameters for sterilization cycles.
- Colour change indicates effective sterilization processes.

200/Pack

Code: (CI1100)

Reg. \$31

Sale: \$25



CI2-P1

Bowie-Dick Test Card


- Bowie Dick Test, Consists of a Type 2, metals free chemical indicator. Used in Class B Sterilizers.
- Use daily to monitor the effectiveness of air removal.
- Cycle information can be recorded on the back of the indicator.

120/Box

Code: (CI1300)

Reg. \$251

Sale: \$201



Bowie-Dick Test Card Holder

- For Bowie Dick Test Card
- Metal Card Holder

Code: (CI1400)

Reg. \$115.50

Sale: \$92.50



Sensibles[™]

COMPATIBLE WITH ALL BRAND SENSORS!!


UNIVERSAL SENSOR HOLDER STARTER KIT

- Everything you need to capture the perfect image is inside!
- Adjusts easily; holds securely
- Perfect for Anterior, Posterior and Bite-Wing images
- Autoclavable and reusable
- Durable and 100% latex-free

Code: (DS1103)

Reg. \$283.92

Sale: \$215⁵⁰



Include:

- Positioning Arm Posterior x1
- Positioning Arm Bite Wing x1, Positioning Arm Anterior x1
- Bite Block Large Size x2, Bite Block Medium Size x2
- Posterior Aiming Ring x1, Anterior/Bite Wing Aiming Ring x1

Using Sensibles is as easy as 1-2-3




1 Attach Positioning Bar



2 Add Aiming Ring



3 Snap In Sensor



For Bite-Wings, adjust vertical position

FREE GOODS REDEMPTION: 1.Go to “<http://beyesdental.com/en/redeems/new>” 2.Fill out the form and submit (Promo code: **CAQ425**)

*Free goods must be of equal or lesser value and cannot be combined with any other offers. All offers will require a proof of purchase dated from October 1 - December 31, 2025 and must be reflected on one invoice. Offers must be redeemed no later than January 15, 2026. Please allow 4 to 8 weeks for promotion to be processed. All brands are holders of their respective trademarks.